WE HAVE REDISCOVERED HEALING DENTISTRY

# LESS *IS* MORE

WWW.UKDENT.COM

LIMITED TIME SPECIAL OFFER ON NEW HEALOZONE x4



#### HealOzone / Fuji Triage: 6 months by Prof Ray Bertolotti



NEW

CAN I AFFORD TO WORK WITHOUT IT?

# NEW EQUIPMENT DECISIONS, DECISIONS

PURCHASE DEDIMONS FIRST TIME, EVERY TIME, FOR EVERY PRACTICE

#### ROI RETURN ON INVESTMENT

NET INCOME (OR BENEFIT) AS A PERCENTAGE OF THE ANYIGINAL AMOUNT INVESTED

FOR EXAMPLE: A CERTIFICATE OF DEPOSIT RETURNING \$50 ON \$1000 INVESTMENT HAS ROI OF 5% PER YEAR Many Dentists report earning the full cost of the HealOzone within one month

Treatment of all cavity preparations with the HealOzone



Specialist in Endodontics, Prosthodontics, Restorative Dentistry and BUPA Consultant in Oral Surgery.

> Head of Dentistry University of Warwick, UK







#### **University of Warwick**



#### Warwick Dentistry



#### **University of Warwick**

- > 20,000 students
- 4,000 from overseas 120 nationalities
- Alumni:125,000 from
- 171 countries
- 290 hectare campus
  - 290 nectare campus





University of Warwick Academic Excellence

Ranked 3<sup>rd</sup> in the University League Table for 2011



"To be a leading medical school in the UK, internationally renowned fo the high quality and relevance of our education programmes and for the excellence and significance of **Warwick Dentistry** 





#### **CPD - Facts and Figures**

- Serve a very wide range of health professionals
- Largest provider in the UK
- >10,000 people have enrolled on one of our CPD programmes
- 1,000 students currently on postgraduate programmes
- Overall numbers are increasing yearly
- Over 20 masters degrees available
- 100 accredited modules

WARWICK

#### Warwick Dentistry

- Orthodontics
- **Lingual Orthodontics**
- Diploma Orthodontic Therapy
- Endodontics
- Prosthodontics
- Periodontology
- Restorative Dentistry
- Oral Surgery
- Implantology
- eMasters in Aesthetic Dentistry Distance
- Learning
  - 3 years full time supervised clinical training

#### One to one mentoring





#### University guide 2011: University league table All UK universities ranked by the Guardian according to teaching excellence

The Guardian University Guide and The Guardian Postgraduate Guide are now available from Guardian Books • Datablog: download the 2011 tables

guardian.co.uk, Tuesday 8 June 2010 00.03 BST

#### University league table More about the university guide Go to:

Rating	Name of institution	Guardian score /100	Satuled overall (%)	Satisfied with teaching (%)	Satisfied with feedback (%)	Student staff ratio	Spend per student (FTE)	Ave entry tariff	Value added score /10	Job after 6 mths
1 (1)	Oxford	100.0				10.8	10	531	6.9	80
2(2)	Cambridge	94,5	92	91	72	11.7	9.59	546	5.8	78
3 (4)	Warwick	81.4	88	87	61	13.3	8.39	463	6.2	75
4 (3)	St Andrews	81.3	92	93	70	13.1	6.02	454	6.9	71





# Janina

Toothpaste Mouthrinse Spray

From Dental Directory or CTS



**Great Oral Health** 

#### Why look for Gum disease ?



#### Why look for Gum disease ?



Feb 28th 2009

#### **Perio-systemic link**



#### Periodontal disease is a risk marker for coronary heart disease

eriodontal disease and coronary heart disease incidence: a systematic review and meta

#### **Conclusion:**

Periodontal disease is a risk factor or marker for CHD that is independent of traditional CHD risk factors, including socioeconomic status.

Further research is warranted in this important area of public health.

#### Periodontal treatment could improve glycaemic control in diabetic patients

Efficacy of periodontal treatment on glycaemic control in diabetic patients: a meta-analysis of interventional studies. *Diabetes Metab 2008; 34:497-506*-**Derre L, Vergnes JN, Courdy P, Stoc** M.

#### Conclusions

The present meta-analysis represents the best information available to date that addresses this issue, and suggests that periodontal treatment could improve glycaemic control.

#### Celebrity Endorsement!

"Smoking kills......

lf you are killed - you`ve lost a

very important part of your

life" Brooke Shields - smoking



 Some examples of minimal invasive dentistry using HealOzone



#### Can we seal in caries?

No 70% of "sealed in" occlusal caries progressed

ACTA Research, Poorterman et al, Caries Research 2003; 37: 29 – 33. Weerheijm KL, de Soet JJ, van Amerongen WE, de Graaff J J Dent Child 1992; 27: 417 – 423

Cariogenic microorganisms found despite sealed restoration with soft and moist dentine, ie active caries under the sealed restorations McDonald SP, Sheiham A.

Caries progressed in 5% of the SDF/  $SnF_2$  group and 11% of the SDF/  $SnF_2$ / mimimal preparation cavities/ composite resin group in 18 months.

Int Dent J 1994, 44, 465 – 470.

#### Erin Mahoney and Michael Swain

No remineralisation beneath GIC or GIC/Ca(OH)<sub>2</sub> or Composite resin "sealants" after 3-4 months

UMIS and SEM-BSE methodology in Archives of Oral Biology 2004







•Some examples of arrested Root Caries follows using Ozone





Selection of relevant references are in the following slides which have been requested Almost 100 publications are accessible on www.iadr.org and the Pan European IADR and divisional web sites, proving the efficacy of the HealOzone







Am J Dent. 2004 Feb;17(1):56-60.

Effect of ozone on the oral microbiota and clinical severity of primary root caries.

Baysan A, Lynch E

Caries Res. 2000 Nov-Dec;34(6):498-501.

Antimicrobial effect of a novel ozone- generating device on microorganisms associated with primary root carious lesions in vitro.

Baysan A, Whiley RA, Lynch E

#### **BENEFITS**

- CONSERVES TOOTH STRUCTURE
- 1 DAY ENDODONTICS

OZON

- CLEANSES PERIO POCKETS EFFECTIVELY with McKenna adaptation
- IMPROVES APICECTOMY SUCCESS
- MORE PREDICATBLE TREATMENT
- FASTER TREATMENT

#### Benefits for the general dentist

- Effective: at least 99.99% of all oral bacteria (incl. Enterococcus Faecalis) are being killed.
- Fast Treatment.
- Easy: Integrable in existing treatment procedure with no need for special measures.
- Safe, with no side-effects: Ozone does not affect the healthy surrounding tissue. No known side-effects

#### **Benefits for the patients**

- Painless
- Root canal treatment in a single visit
- Improves implant augmentation success rate
  "Inflammation free" dental pocket after perio
- Treatment

   • Less risk of pulpal exposure in treating deep decay
- Less risk of pulpar exposure in treating deep decay
   Less trauma for children and dental phobics
- Reduces the need for local and systemic antimicrobials



Benefits: • Shortens curing time by over 80% • Improves composite flow by 68% • Reduces Microleakage • Improves physical properties









if you are going to watch something you have to be very sure of what you are watching









#### The CarieScan PRO



#### How the PRO works

- The product platform is based on the application of a technique called ac impedance spectroscopy (ACIST).
- The measurement relies on the application of small electrical signals through the tooth while monitoring the response to the sensor.
- Analysis of the impedance response relates to the physical structure of the tooth



# World Beating Accuracy

94.8% accurate at detecting caries

- Fewer than 7% false positives and 7 times less than our nearest competitor
- Increase your and your patients confidence that the correct diagnosis has been made



-

#### Simple to use.....

Simply air dry the tooth as you would for a visual examination

- · Place the sensor on the tooth/site
- A measurement will be displayed immediately
- To take another reading simply place the sensor on another tooth/site
- No calibration required

#### ....and understand

- The PRO has a simple 0 -100 scale
- · 0-50 Low probability of caries
- 51-90 Medium probability of caries
- 91-99 High probability of caries
- 100 Operative Care Advised
- Quantifiable and repeatable output allows for on-going patient monitoring and an increase in repeat preventative treatments





Д

105

Preparation work can be passed to Hygienists saving time and money

More preventive and restorative work

#### **Clinicians Report**



- 'Easiest and fastest to use'
- 'Performed better than all previous products in rigorous trials'
- · The only product to be marked as 'excellent' in any categories (2)





#### **Testimonials**

What sets the CarieScan PRO apart from other units is its improved sensitivity and specificity with an unmatched accuracy. This high degree of accuracy means fewer false positives, giving me more confidence in my diagnosis'

Dr MICHAEL MIYASAKI, DDS (Sacramento, CA, USA)

"I can confidently tell you that among the 11 devices we have in our caries detection clinic, CarieScan is the most reliable and accurate device for detection and monitoring of initial caries lesion on all surfaces."

PROF BENNETT TAMAECHI BDS, MS (SATX, USA)

#### **Testimonials**

'I use a Cariescan PRO daily. I have found it to be an invaluable device. The PRO has been 100% reliable for the detection of early carious lesions and is now routinely used in place of radiographs' IAN ROBERTSON BDS (UK)

'It is so simple, patients have responded well to it and it has increased visits to the dentists in the practice' MHARI COXON, RDH (UK)

#### Edward Lynch

Evidence-based efficacy of ozone for root canal irrigation.

J Esthet Restor Dent. Oct; 2008;20(5):287-293.



Efficacy of ozone on survival and permeability of oral microorganisms

Nagayoshi M, Fukuizumi T, Kitamura C, Yano J, Terashita M, Nishihara T.

Oral Microbiology and Immunology Volume 19 Page 240-246 - 2004.

Antimicrobial Effect of Ozone on Bacteria Invading Dentinal Tubules

Nagayoshi M., Kitamura C., Fukuizumi T., Nishihara T., Terashita M.

Journal of Endodontics, Vol. 30, No.11, 2004, 778-781

# Mechanism of Action

<u>J Dent Res.</u> 2002 Jun;81(6):422-7. 1H and (13)C NMR spectroscopic analysis of human saliva.

Silwood CJ Lynch E Claxson AW Grootveld MC

J Biol Inorg Chem. 2002 Jan;7(1-2):46-57.

1H NMR investigations of the molecular nature of lowmolecular-mass calcium ions in biofluids.

Silwood CL, Grootveld M, Lynch E.

#### Molecular mechanisms

Silwood C.J.L., Lynch E., Seddon S., Sheerin A., Claxson A.W.D. And Grootveld M. (1999)

<sup>1</sup>H NMR Analysis of Microbial-Derived Organic Acids in Primary Root Carious Lesions and Saliva.

NMR in Biomedicine 12: 345-356.

#### <sup>1</sup>H NMR spectra









The use of ozone and 0.2% chlorhexidine in the treatment of periodontitis patients: A clinical and microbiologic study

Kshitish D, Laxman VK. Indian J Dent Res. 2010 Jul-Sep;21(3):341-8. Effect of ozone on periodontopathogenic species

S Eick, M Tigan, and A Sculean

IADR SAN DIEGO MARCH 2011, PAPER 3070 Results showed a high efficacy of ozone against Fusobacterium nucleatum, Porphyromonas gingivalis and Aggregatibacter actinomycetemcomitans.

Most of the strains in a concentration of 10<sup>5</sup> were completely eliminated after two-fold 18 s application of ozone Effectiveness of ozone against periodontal pathogenic microorganisms

Huth et al, in press 2012

Compared to CHX 0.2%, aqueous ozone 20 µg/ml and gaseous ozone down to 4 g/m3 showed a significantly higher effectiveness.



Is Ozone useful for Periimplantitis management or prevention?

#### Dr Dan Mc Kenna

- Proven successful prevention of Peri-implant mucositis using Ozone
- Awarded First Prize IADR 2009

MPhil Thesis 2010

Peri-implant mucositis is a significant clinical problem. Huth et al (J Dent Res. 2007;86 (5):451-456) reported potential benefits of using Ozone in periodontal treatment. As ozone is the most powerful antimicrobial agent we could use in dentistry it seemed appropriate to assess if Ozone could prevent peri-implant mucositis (McKenna et al PEF 2008).

The aim of this study was to assess if metrology could quantify experimentally produced peri-implant mucositis with and without HealOzone treatment.

The Local Research Ethics Committee (LREC) granted ethical approval and informed consent was also given by all patients after reading the study information sheet.

#### Methods:

Twenty subjects were enrolled in a randomized, double – blind, controlled, single center study after a pre-trial phase to achieve clinically healthy gingivae

All subjects at the end of the pretrial phase were asked to refrain from tooth brushing that area by wearing a gum shield over the chosen implants whilst tooth brushing.

Each of these patients had two implants which were randomly allocated as either test or control. Subgingival administration of ozone (HealOzone, KaVo) and saline (test implant site) or air and saline (control implant site) was delivered subgingivally for 60 seconds on day 0, 7 and 14. Plaque, bleeding and modified gingival indices and impressions were recorded by a single operator at 0 and 21 days in this 21 day experimentally induced periimplant mucositis study

The patients and dentist (investigator) were unaware of the treatments applied to each implant site

Special additional switch to deliver air instead of ozone (Rear view)





HealOzone handpiece adapted by Mc Kenna to treat subgingivally with HealOzone silicone cup and cannula

Mode of HealOzone handpiece cap applications





Paper point to judge depth of gingival crevice



Stacking of silicone cups



3D Scanning software showing crosssectional slice of stone model



Results: There were significant differences in plaque (p-value <0.01), gingival (p-value = 0.011) and bleeding (p-value <0.01) index scores between the treatments with ozone and saline showing the optimal tissue health scores.

Plot of the mean plaque index for each treatment







The bleeding index at 21 days was less for the ozone treated sites than for the air treated sites (p<0.01) and the change in the gingival volume from the laser scanned impressions correlated with this result.







# Can HealOzone whiten teeth?



<u>J Ir Dent Assoc.</u> 1994;40(3):78-82. Safety aspects concerning the therapeutic and cosmetic applications of hydrogen peroxide (H2O2)--containing gels, whiteners, oral rinses and dentifrices. Lynch E, et al

<u>J Esthet Dent.</u> 1991 Sep-Oct;3(5):162-8. Is home tooth bleaching gel cytotoxic? <u>Tse CS, Lynch E, Blake DR, Williams DM.</u> Clin Oral Investig. 2010 Feb;14(1):1-10. Undesirable and adverse effects of toothwhitening products: a review.

Goldberg M, Grootveld M, Lynch E.

Genotoxicity and carcinogenicity only occur at concentrations that are never reached during dental treatments.

# Mechanism of Action

<u>J Ir Dent Assoc.</u> 1995;41(4):94-102.

Molecular mechanisms of the bleaching actions associated with commercially-available whitening oral health care products.

Lynch E, et al.

#### Mechanism of Bleaching

Oxidation

Especially from free radical generation with superoxide (O·) radicals and hydroxyl radicals (OH·) J Dent Res. 2002 Jun;81(6):422-7.

1H and (13)C NMR spectroscopic analysis of human saliva.

Silwood CJ Lynch E Claxson AW Grootveld MC

J Biol Inorg Chem. 2002 Jan;7(1-2):46-57.

1H NMR investigations of the molecular nature of lowmolecular-mass calcium ions in biofluids.

Silwood CL, Grootveld M, Lynch E.

#### Molecular mechanisms

Silwood C.J.L., Lynch E., Seddon S., Sheerin A., Claxson A.W.D. And Grootveld M. (1999)

<sup>1</sup>H NMR Analysis of Microbial-Derived Organic Acids in Primary Root Carious Lesions and Saliva.

NMR in Biomedicine 12: 345-356.

# **Activators**

Oper Dent. 2011 Mar-Apr;36(2):162-8. **The influence of chemical activation on tooth bleaching using 10% carbamide peroxide**. <u>Batista GR, Barcellos DC, Torres CR, Goto EH,</u> <u>Pucci CR, Borges AB</u>.

 Adding 0.01% manganese gluconate to 10% carbamide peroxide bleaching gel increased the degree of tooth bleaching after a sevenday treatment and did not influence the resulting shade after 14 days.

Dental Materials Journal 2011; 30(5): 723-729

Effect of light units on tooth bleaching with visible-light activating titanium dioxide photocatalyst

Ayaka KISHI<sup>1</sup>, Masayuki OTSUKI<sup>1</sup>, Alireza SADR<sup>2</sup>, Masaomi IKEDA<sup>3</sup> and Junji TAGAMI<sup>12</sup>

This study evaluated the influence of different light sources on the efficiency of an office bleaching agent containing visible-light activating titanium dioxide photocatalyst (VL-TiO2). Dental Materials Journal 2011; 30(5): 723–729 Effect of light units on tooth bleaching with visiblelight activating titanium dioxide photocatalyst Ayaka KISHI et al

The blue and violet LED light curing was more effective than blue LED light curing and conventional halogen lamp for tooth bleaching with an office bleaching agent containing visible-light activating titanium dioxide photocatalyst (VL-TiO2)

#### **Chlorine dioxide**

Lynch E. *et al.* (1997)

Multicomponent spectroscopic investigations of salivary antioxidant consumption by an oral preparation containing the stable free radical species chlorine dioxide (CIO<sub>2</sub>\*).

Free Radical Research 26:209-237.

#### Peroxoborate

Lynch E . *et al.* (1999)

Multicomponent evaluations of the oxidising actions and status of a peroxoborate-containing tooth-whitening system in whole human saliva using high resolution proton NMR spectroscopy.

*Journal of Inorganic Biochemistry* 73: 65-84.

#### **Surgery Bleaching**

Ozone activates hydrogen peroxide Can HealOzone whiten teeth?

**Ozonated gel and Composite restorations** 

n Dental Journal 2010; 55: 390-398

dor: 10.11116 1654-7819 2010

SCIENTIFIC ARTICLE

Effect of bleaching versus repolishing on colour and surface topography of stained resin composite

M Abd Elhamid,\* R Mosallam\*

\*El Seka El Haded Hospital, Cairo, Egypt.
†Department of Operative Dentisity, Faculty of Oral and Dental Molicine, Cairo University, Egypt.

#### **Ozonated gel and Composite Restorations**

"Superior whitening effect was demonstrated with the ozonated gel. Ozonated gel showed statistically significant lowest roughness compared to both carbamide peroxide and polishing paste."

#### **Ozonated gel and composite restorations**

Whitening nethod (I)	Whitening method (J)	Mean difference (I-J)	Р
Carbamide peroxide	Polishing paste	-0.800	0.208
	Ozonated gel	*1.600	0.003*
Polishing paste	Carbamide peroxide	0.800	0.208
01	Ozonated gel	*2.400	0.000*
Ozonated gel	Carbamide peroxide	*-1.600	0.003*
	Polishing paste	*-2,400	0.000*

Ozone is an efficient bleaching agent with the least

adverse effect on surface roughness."





Bleaching; Ozone and hydrogen peroxide combinations Effect of ozone, chlorine and hydrogen peroxide on the elimination of colour in treated textile wastewater by MBR.

> Water Sci Technol 2004;49(4):299-303

Brik M; Chamam B; Schoberl P; Braun R; Fuchs W Effect of oxygen, ozone and hydrogen peroxide bleaching stages on the contents and composition of extractives of Eucalyptus globulus kraft pulps. Bioresour Technol 2006 Feb;97(3):420-8

Freire CS; Silvestre AJ; Pascoal Neto C; Evtuguin DV

Dechlorination of chlorophenols found in pulp bleach plant E-1 effluents by advanced oxidation processes.

> Bioresour Technol 2005 May;96(8):897-906

Wang R; Chen CL; Gratzl JS

Rate of dibutyIsulfide decomposition by ozonation and the O3/H2O2 advanced oxidation process.

J Hazard Mater 2009 May 30;164(2-3):1364-71

Popiel S; Nalepa T; Dzierzak D; Stankiewicz R; Witkiewicz Z Dehalogenation, degradation and mineralization of diuron by peroxone (peroxide/ozone) treatment.

J Environ Sci Health A Tox Hazard Subst Environ Eng 44;6(630-8 S1093-4529

Catalkaya EC; Kargi F

Kinetics study on photochemical oxidation of polyacrylamide by ozone combined with hydrogen peroxide and ultraviolet radiation. J Environ Sci (China)

2006;18(4):660-4

Ren GM; Sun DZ; Chung JS

Elimination of Listeria monocytogenes biofilms by ozone, chlorine, and hydrogen peroxide.

J Food Prot 2005 Mar;68(3):494-8

Robbins JB; Fisher CW; Moltz AG; Martin SE An ozone/hydrogen peroxide/microwave-enhanced advanced oxidation process for sewage sludge treatment.

J Environ Sci Health A Tox Hazard Subst Environ Eng 42;8(1177-81 S1093-4529

Yin G; Liao PH; Lo KV

Oxidative degradation of Nnitrosodimethylamine by conventional ozonation and the advanced oxidation process ozone/hydrogen peroxide.

Water Res 2007 Feb;41(3):581-90

Lee C; Yoon J; Von Gunten U

Treatment of volatile organic chemicals on the EPA Contaminant Candidate List using ozonation and the O3/H2O2 advanced oxidation process. Environ Sci Technol 2006 Apr 15;40(8):2734-9

Chen WR; Sharpless CM; Linden KG; Suffet IH Grootveld M, Silwood C, Lynch E.

Biofactors 2006; 27: 5 – 18.

Oxidative Consumption of Biomolecules by Ozone; Clinical Relevence.

J Ir Dent Assoc. 1996;42(3):74-6.

Professional bleaching of teeth in dental practice techniques.

Samarawickrama DY, Lynch E, et al

### **Surgery Bleaching**

Ozone activates hydrogen peroxide







## Zero Sensitivity Perfect Tray (Inventor, Dr. Wyman Chan, UK Patent No. 2 416 310, 2 445

- 298)
- Non-Reservoir (CRA 1997 Gordon Christensen)
- Non-Scalloped, extend 2-3 mm over gingival . margins, preventing saliva entering the tray to dilute the gel (Haywood 2002)



Perfect Seals, like post dams, on the inside labial preventing the gingival crevicular fluid (GCF) and saliva contaminating and consuming the bleaching agent and also sealing the gel within the tray







- The dimples serve two purposes: The location marks and
   the amount of bleaching gel to be loaded inside the tray
- Only a small amount of gel is required



Bleaching gel is prevented from getting in contact with exposed dentine and oral soft tissues by the Perfect Seals



Bleach only enamel and gel is locked in by the Perfect seal (enamel is the only non-vital tissue in our body when fully matured) Oxidants (gas) are pushed back into the enamel pores enhancing the bleaching effect



- Minimal amount of gel needed
- 75mg (6%) per upper application (502mg) Zero sensitivity on teeth and oral soft •
- tissues.



#### Smile Studio Perfect Tray HB

- Carbamide Peroxide (CP) 10%-16% -Overnight wear
- CP 18% 1-2 hours twice a day or overnight
- $H_2O_2$  (HP)  $\frac{1}{2}$  hour twice a day
- 2 to 6 weeks
- Lower wear time if experience sensitivity







- Gentle raised temperature of ō bleaching gel at gel/enamel interface
- Enhancing diffusion of gel inside enamel wy10 10% H<sub>2</sub>O<sub>2</sub>









### **Internal bleaching**

HealOzone and the walking bleach technique

#### Access and Seal over GP

37% Phosphoric Acid Etch

Paste of Sodium Perborate Powder and  $H_2O_2$ 

HealOzone 30 seconds

Seal with GIC 1 wk

### Biocompatible











healOzone X4Curozone

•www.ukdent.com •E.lynch@warwick.ac.uk

HealOzone handpiece from CurOzone Germany









#### Treatment of dentine with HealOzone

HealOzone





Knight GM et al., Aust Dent J 2008;53:349-353.

# Ozone (O<sub>3</sub>) treatment of caries:

- The "ecological niche" of acidogenic and aciduric microorganisms is eliminated, for at least 14 weeks
- Remineralization overpowers
   demineralization
- Rapid remineralization occurs

•Some examples of radiograhic changes in Ozone treated caries follows













"I have included a sequence of Xrays I did on one of these patients.'

X-ray 1: Child patient aged 4 years. Presents with rampant bottle decay on all teeth. No pain or infection present. Suggested GA with pulpotomy to try to save 55 to allow 16 to erupt into favourable position. Parents refused general . anesthesia.

X-ray 2: Child presents 6 months later. Again, no pain or infection. Parents do not want GA. Child uncooperative but allowed ART technique with Fuji IX.





X-ray 3: Child returned over one year later as an emergency with toothache. ART technique had failed on 55. Child still unable to have conventional dentistry performed in the chair. Parents still unwilling to accept GA. Ozone from the Healozone and Fuji VII was used.

X-ray 4: Child came for a "check-up" just over 1 year later. Miraculous! No pain or infection. The 16 has erupted perfectly into position.

My question here is: how is the tooth "healing"? Please can you explain the mechanism by which this hyper-mineralised tissue is formed?

Huth K, Paschos E, Brand and Hickel R. Effect of Healozone on non cavitated fissure carious lesions in permanent molars. A controlled prospective study.

**American Journal of Dentistry** 2005; 18: 223 - 270.

Most Dentists charge the same fee for Ozone treatment and sealing (flowable composite) of an occlusal caries lesion as they do for a posterior composite. This saves the Dentist time

Reversal of deciduous caries after HealOzone treatment Dahrhardt J, Jaeggi T, Lussi A. Treating open carious lesions in anxious children with ozone. A prospective controlled clinical study. American Journal of Dentistry 2006; 19: 267 – 270.

Reversal of Deciduous Caries using HealOzone

> M Phil Thesis, UK 2004

> > **OT.Abu-Salem**

HealOzone uses with Caries

Professor Sebastian Ciancio, USA Biological Therapies in Dentistry 2005 "Studies with over 2000 Patients have shown remineralisation of Healozone treated permanent and deciduous occlusal surfaces as well as root caries". "The HealOzone is a novel new and painless way to treat early caries".

Amna Al Shamsi

- PhD thesis 2007
- HealOzone significantly reduces caries incidence around orthodontic brackets (9% versus 28%)

Dahrhardt J, Jaeggi T, Lussi A. Treating open carious lesions in anxious children with Healozone. A prospective controlled clinical study. American Journal of Dentistry 2006; 19: 267 – 270.

Huth K, Paschos E, Brand and Hickel R. Effect of Healozone on non cavitated fissure carious lesions in permanent molars. A controlled prospective study.

American Journal of Dentistry 2005; 18: 223 – 270. Effect of Healozone treatment on different cariogenic microorganisms in vitro.

Swed Dent J 2008;32(3):139-47 (ISSN: 0347-9994) Fagrell TG; Dietz W; Lingstrom P; Steiniger F; Noren JG The Use of Ozone in Dentistry and Medicine. Part 1.

Baysan A and Lynch E.

Primary Dental Care, 12; 2: April 2005, 47 – 52.

The use of Ozone in Dentistry and Medicine.

• Part 2.

- Baysan and Lynch
- Primary Dental Care 2006; 13: 37 – 41

Anti-microbial effects of a novel ozone generating device on micro-organisms associated with root carious lesions

A. Baysan, R. Whiley and E. Lynch Caries Research 2000;34:498-501.

99% microbial killing achieved after ozone treatment HealOzone treatment for deep caries or as an alternative to stepwise excavation

Studies in London and Isle of Wight Reversal of root caries using HealOzone - A 12 month longitudinal study

A. Baysan (London)

PhD Thesis University of London 2003

#### Am J Dent 2004; 17: 54 – 60.

- Baysan A and Lynch E.
- 91% Reversal of Caries at 5.5 months after Healozone treatment.
- 99% Reduction in Microorganisms at 5.5 months.

•The Root Caries deeper then 2mms did not reverse. Gerodontology 2003; 20: 106 – 114.

Clinical reversal of root caries using Healozone, double blind, randomised, controlled 18-month trial.

- Julian Holmes
- 100% Reversal of Non Cavitated Root Caries at 18 months

Antibacterial effect of Healozone on cariogenic bacterial species.

J Dent. 2009 Jun;37(6):449-53

Johansson E Claesson R van Dijken JW

The influence of Healozone on microleakage and fissure penetration of different sealing materials.

Coll Antropol. 2009 Mar;33(1):157-62.

Dukić W, Dukić OL, Milardović S.

Ozone improves lipopolysaccharideinduced responses of an odontoblast-like cell line.

J Endod. 2009 May;35(5):668-72

<u>Noguchi F, Kitamura C, Nagayoshi M,</u> <u>Chen KK, Terashita M, Nishihara T</u>.

The inability of Streptococcus mutans and Lactobacillus acidophilus to form a biofilm in vitro on dentine pretreated with Healozone.

Aust Dent J. 2008 Dec;53(4):349-53

Knight GM, McIntyre JM, Craig GG, Mulyani, Zilm PS. Treating sensitive cervical areas with Healozone. A prospective controlled clinical trial.

Am J Dent. 2008 Apr;21(2):74-6

<u>Dähnhardt JE, Gygax M, Martignoni B,</u> <u>Suter P, Lussi A</u>. Skaug N, Strand G, Nielsen O

**ORCA & Caries Research 2006** 

85% Reduction in Mutans Streptococci with HealOzone and 78.6% Reduction in Lactobacilli

#### <u>Khairul Matin</u>, Junji Tagami

99% Microbial killing with HealOzone

Dept of Restorative Sciences and Cariology Tokyo Medical and Dental University Professor Michael Noack and Suzanne Kneist

**ORCA & Caries Research 2006** 

99% reduction in microorganisms following 20 seconds of HealOzone treatment HealOzone treatment for deep caries or as an alternative to stepwise excavation

Studies in London and Isle of Wight

#### Are Ozone systems safe?

Ozone air levels from a dental ozone gas delivery system

Johansson E, Andersson-Wenckert I, Hagenbjork-Gustafsson A, van Dijken JWV

Acta Odontol Scand. 2007 Nov;65(6):324-330

#### Conclusion

#### HealOzone is safe

# Assessment of the safety of two ozone delivery devices

Millar BJ, Hodson N.

J Dent. 2007 Mar;35(3):195-200

#### Conclusion

HealOzone is safe and the Ozone system which blows out Ozone is not safe J Dent Res 86 (5): 451–456, 2007 Huth KC et al. Effect of Aqueous Ozone on the NF-kB System

Ozone exerts inhibitory effects on the NF-kB system suggesting it has an antiinflammatory capacity. IkBa proteolysis, cytokine expression and kB-dependent transcription were prevented.

#### Huth et al.

Effect of ozone on oral cells compared to established antimicrobials.

Eur J Oral Sci 2006; 114: 435 - 440.

Filippi A

The effects of ozone on epithelial wound healing.

Dtsch Zahnarztl Z 2001; 56: 104 – 108.

Research Awards to Ozone researchers

• Dr Dan Mc Kenna	• Dr Jameela Alawadi	Professor Martin Grootveld
<ul> <li>Proven successful prevention of Peri-implantitis using HealOzone</li> <li>Awarded First Prize IADR</li> </ul>	<ul> <li>Proven successful use of HealOzone in root canal therapy</li> </ul>	<ul> <li>Proven successful management of Dental Unit Water Lines using Ozone</li> </ul>
2009 (International Association for Dental Research annual meeting)	Awarded First Prize IADR     2008	• Awarded First Prize IADR 2007

Professor Dr Liviu Steier	• Dr Wyman Chan	• Dr Julian Holmes
<ul> <li>Proven successful use of HealOzone and sealing to manage caries</li> </ul>	<ul> <li>Awarded First Prize IADR for his HealOzone research 2005</li> </ul>	• Awarded First Prize IADR for his HealOzone research 2004
Awarded First Prize IADR     2006		

• Dr Aylin Baysan	• Dr Layla Abu-Naba'a	• Helene Domingo
• Awarded two First Prizes at IADR	• Awarded the prestigious Basil Bibby cariology award at IADR partly for her HealOzone research 2001	<ul> <li>Awarded First Prize IADR for HealOzone research</li> </ul>
- Amna Al Shamsi
- PhD thesis 2007
- HealOzone significantly reduces caries incidence around orthodontic brackets (9% versus 28%)
- Ozone therapy in the treatment of avascular bisphosphonate – related jaw osteonecrosis
- Agrillo et al. J Craniofac Surg 2007; 18: 1071 – 1075
- Ozone increased the complete healing of lesions with the disappearance of symptoms and brings lesion progression down to zero. Ozone's "benefits were

Objective Quantitative Use of Saliva to Reflect Health or Disease

> E Lynch IADR Symposium 2009





Disease	Device	Stage	Funding
Congestive Heart Failure	Blue Scale	Clinical trial complete	۲
Aortic Aneurysm	Blue Box + Sphygmomanometer	Design completed	DeBakey Center
Asthma	Single Breath Spiromter + mobile phone+video game + social network	Under construction	Challenge Grant (applying)
Diabetes/Metaboli c Syndrome	Tunable laser spectrometer & hand held breath detector	Basic science development	Metholist
HIV/AIDS	Medication dispenser/mobile phone/social network	Planned/awaitin g funding	GATES (medicine) (applied)
Buruli Ulcer	Mobile phone camera	Software testing	

### ADVANTAGES OF SALIVA AS BIOFLUID MEDIUM FOR DIAGNOSTIC PURPOSES

- Ease of Collection
- Low cost of Collection
- Non-invasive Collection for Patients (reduces *Anxiety* and *Stress*)
- Facilitates Collection of Multiple Samples for Time-Dependent Monitoring Purposes
- Easy to Handle and Deal with

Markers for pregnancyrelated disorders, such fetal aneuploidy, preterm birth, preeclampsia, intraamniotic infection and fetal stress. MANY OTHER DIAGNOSTIC ANALYTES ('BIOMARKERS') HAVE BEEN SHOWN TO BE PRESENT IN HUMAN SALIVA

For example:

- Steroid Hormones
- HIV Antibody
- Those for Hepatitis A, B and C

Unique diagnostic panels of salivary mRNAs in subjects with Sjögren's disease.

Four salivary mRNAs (OAZ, SAT, IL8, and IL1b) collectively have a discriminatory power of 91% sensitivity and specificity for oral cancer detection.

### **Recommended Reading**

The Axelsson Series on Preventive Dentistry

Ozone – the Revolution in Dentistry. Quintessence Dec 2004 New opinions are always suspected, and usually opposed, without any other reason but because they are not already common

John Locke 1690





E.LYNCH@WARWICK.AC.UK



E.LYNCH@WARWICK.AC.UK

See

www.realityesthetics.com

**EL103** 

E.LYNCH@WARWICK.AC.UK



E.LYNCH@ WARWICK.AC.UK **Ozone:** he Revolution in Dentistry







### E.LYNCH@WARWICK.AC.UK

THE UNIVERSITY OF



- Effect of Ozone on Dental Caries
   Progression
- Seghi et al
- AADR Dallas April 2008
- Ozone significantly (P=0.003) reduced caries progression in experimental rats

Oxidation of Biomolecules by Ozone

- E Lynch et al
- AADR Dallas April 2008
- Oxidation of Cysteine and Methionine proving Ozone (TherOzone) can combat oral malodour

- Bleaching of compounds responsible for Tooth Discoloration by Ozone
- H Domingo et al
- AADR Dallas 2008
- Ozonated water (TherOzone) bleaches

### Han S. Uhm, Kwang H. Lee, and Baik L. Seong.

"Inactivation of H1N1 viruses exposed to acidic ozone <u>water</u>."

*Applied Physics Letters* 95, 173704 (2009).

Ozone use in RCT, benefits of Ozone for healing, periodontology and safety.

### HealOzone is safe

### **Ozi-cure is not safe**

Miller B and Hodson N, Assessment of the safety of two ozone delivery devices. J Dent 2007; 35: 195 – 200. Ozonated water improves lipopolysaccharide-induced responses of an odontoblast-like cell line.

J Endod. 2009 May;35(5):668-72

<u>Noguchi F Kitamura C Nagayoshi M</u> <u>Chen KK Terashita M Nishihara T</u> Ozone air levels from a dental ozone gas delivery system

Acta Odontol Scand Nov 2007.

Johansson E, Andersson-Wenckert I, Hagenbjork-Gustafsson A, van Dijken JWV.

HealOzone is safe.

<ul> <li>Wounds receive more Oxygen when Ozonated Water is applied.</li> <li>Enhances the phagocitary activity defending cells</li> </ul>	Ozone therapy in medicine and dentistry.	<ul> <li>Scientific and medical aspects of Ozone therapy. State of the art.</li> </ul>
<ul> <li>Accelerates migration of epithelial cells</li> <li>Activates Fibroblasts</li> </ul>	J Contemp Dent Pract 2008;9(4):75-84 (ISSN: 1526-	Professor Velio Bocci
<ul> <li>Increases collagen synthesis</li> <li>Improves cell proliferation</li> <li>Increases chemotaxis of monocytes and</li> </ul>	3711)	• Arch Med Res 2006; 37: 425 – 435.
fibroblasts <ul> <li>Increases synthesis of extracellular matrix</li> </ul>	Nogales CG; Ferrari PH; Kantorovich EO; Lage-Marques	<ul> <li>Beneficial effects of Ozone discussed</li> </ul>
Professor Bocci / Professor Filippi	JL	

<ul> <li>Ozone accelerates wound healing</li> </ul>	<ul> <li>Ozone activates cellular metabolism. Masui 1993; 42: 2 – 6.</li> </ul>	J Dent Res 86 (5): 451–456, 2007 Huth KC et al. Effect of Aqueous Ozone on
• Chronic leg ulcers. Aust Farm Physician 1985; 14: 292 - 298	<ul> <li>Ozone raises intracellular ATP. Toxicology 1991; 70: 195 – 202.</li> <li>Ozone increases cytokines relevant to wound bealing especially.</li> </ul>	the NF-kB System Aqueous Ozone exerts inhibitory effects on the NF-kB system suggesting it has an anti-inflammatory capacity.
• Skin grafts. Langenbecks Arch Chir 1995; 42: 2 – 6.	Transforming Growth Factor. Professor Bocci et al, J Biol Regul Homeost Agents 1994; 8: 108 – 112.	IkBa proteolysis, cytokine expression and kB-dependent transcription were prevented.

Evidence-based efficacy of ozone for root canal irrigation.

Lynch E.

J Esthet Restor Dent. 2008;20(5):287-93.

Efficacy of calcium hydroxide, Er:YAG laser or gaseous ozone against Enterococcus faecalis in root canals.

Am J Dent. 2009 Feb;22(1):14-8.

<u>Noetzel J Nonhoff J Bitter K Wagner</u> J <u>Neumann K Kielbassa AM</u> Bactericidal effect of KTP laser irradiation against Enterococcus faecalis compared with gaseous ozone: an ex vivo study.

Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 May;107(5):e73-9

Kuştarci A, Sümer Z, Altunbaş D, Koşum <u>S</u>.

<ul> <li>Jameela Mohammed Alawadi</li> </ul>	• Experience in ozone use for root canal therapy	<u>Nagayoshi M, Kitamura C,</u> Fukuizumi T, Nishihara T,
PhD thesis 2008	• Berukova et al	<u>Terashita M.</u>
<ul> <li>Successful use of Ozone in root canal therapy</li> </ul>	• Stomatologia 2005; 84: 20 – 22.	Antimicrobial effect of ozonised water on bacteria invading dentinal tubules.
• World First Prize IADR 2008	• Proved efficacy of Ozone for root canal therapy.	J Endod 2004; 30: 778 – 781.

Filippi A	<ul> <li>The influence of Ozonated Water on the epithelial wound healing process in the oral cavity</li> </ul>	Huth et al.
The effects of ozonised water on epithelial wound healing.	<ul> <li>Professor A Filippi</li> <li>Use of Ozonater Water clearly</li> </ul>	Effect of ozone on oral cells compared to established antimicrobials.
Dtsch Zahnarztl Z 2001; 56: 104 – 108.	showed an acceleration of wound healing within the first 48 hours, resulting in earlier epithelial would closure after 7 days.	Eur J Oral Sci 2006; 114: 435 – 440.

<ul> <li>Braz Dent J. 2006;17(2):134-8.</li> <li>Antimicrobial potential of ozone in an ultrasonic</li> </ul>	Ozone sterilises 10 <sup>6</sup> cfu Enterococcus Faecalis	Professor Beer and Liviu Steier
cleaning system against Staphylococcus aureus.	Chang H	Witten University
<u>Estrela C</u> <u>Estrela CR</u> <u>Decurcio Dde A</u> <u>Silva</u> <u>JA</u> <u>Bammann LL</u>	IADR 2003	Sterilisation of root canals using HealOzone 2006

Efficacy of Ozone on Survival and Permeability of Oral Microorganisms

M Nagayoshi, T Fukuizumi, C Kitamura, J Yano, M Terashita, T Nishihara

Oral Microbiology and Immunology 2004: 19: 240 – 246. Therapeutic effects of topical application of ozone on acute cutaneous wound healing.

J Korean Med Sci. 2009 Jun;24(3):368-74. Epub 2009 Jun 12

 Kim HS
 Noh SU
 Han YW
 Kim KM
 Kang

 H
 Kim HO
 Park YM

Superficially, longer, intermittent ozone therapy in the treatment of the chronic, infected wounds.

Ortop Traumatol Rehabil. 2003 Oct 30;5(5):652-8.

Białoszewski D Kowalewski M

Dent Today. 2010 Feb;29(2):130, 132-3. Leave decay in my cavity? You must be kidding!

Knight GM, McIntyre JM, Craig GG, Mulyani. Clin Oral Investig. 2010 Jan 7. [Epub ahead of print] Influence of ozone on the composite-to-composite bond.

<u>Magni E, Ferrari M, Papacchini F,</u> <u>Hickel R, Ilie N</u>. J Adhes Dent. 2009 Aug;11(4):287-92. Enamel and dentin bond strength following gaseous ozone application.

<u>Cadenaro M, Delise C, Antoniollo</u> <u>F, Navarra OC, Di Lenarda R,</u> <u>Breschi L</u>.

Monogr Oral Sci. 2009;21:156-63. Epub 2009 Jun 3. Novel preventive treatment options.

Longbottom C, Ekstrand K, Zero D, Kambara M. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 May;107(5):e73-9. Bactericidal effect of KTP laser irradiation against Enterococcus faecalis compared with gaseous ozone: an ex vivo study.

Kuştarci A, Sümer Z, Altunbaş D, Koşum S. J Endod. 2009 May;35(5):668-72. Ozonated water improves lipopolysaccharide-induced responses of an odontoblast-like cell line.

> <u>Noguchi F, Kitamura C,</u> <u>Nagayoshi M, Chen KK,</u> <u>Terashita M, Nishihara T.</u>

Coll Antropol. 2009 Mar;33(1):157-62. The influence of Healozone on microleakage and fissure penetration of different sealing materials.

Dukić W, Dukić OL, Milardović S.

J Dent. 2009 Jun;37(6):449-53. Epub 2009 Apr 1. Antibacterial effect of ozone on cariogenic bacterial species.

Johansson E, Claesson R, van Dijken JW. Am J Dent. 2009 Feb;22(1):14-8. Efficacy of calcium hydroxide, Er:YAG laser or gaseous ozone against Enterococcus faecalis in root canals.

<u>Noetzel J, Nonhoff J, Bitter K,</u> <u>Wagner J, Neumann K,</u> <u>Kielbassa AM</u>.

Stomatologiia (Mosk). 2008;87(6):24-6. [Application of medical ozone in endodontic practice] [Article in Russian]

Bezrukova IV, Petrukhina NB, Dmitrieva NA, Snegirev MV. Some published research with the TherOzone

Some published research with the TherOzone

New therapeutic strategies for the treatment of difficult wounds

G Chir. 2008 May;29(5):212-20

<u>Onesti MG</u> <u>Bitonti A</u> <u>Fino P</u> <u>Ciotti M</u> <u>Scuderi N</u> Beneficial effects of pro-/antioxidant-based nutraceuticals in the skin rejuvenation techniques.

Cell Mol Biol (Noisy-le-grand). 2007 Apr 15;53(1):94-101.

<u>de Luca C</u> <u>Deeva I</u> <u>Mikhal'Chik E</u> <u>Korkina L</u> Therapeutic effects of topical application of ozone on acute cutaneous wound healing.

J Korean Med Sci. 2009 Jun;24(3):368-74

Kim HS, Noh SU, Han YW, Kim KM, Kang H, Kim HO, Park YM.

[Experimental-morphological study of the anti-inflammatory action of ozone-perfluorane complex application]

Stomatologiia (Mosk). 2008;87(2):4-9 Grigor'ian AS, Grigor'iants LA, Guchetl' MN. The case for oxygenozonetherapy.

<u>Bocci V</u>.

Br J Biomed Sci. 2007;64(1):44-9

Periradicular repair after two-visit endodontic treatment using two different intracanal medications compared to single-visit endodontic treatment.

Braz Dent J. 2007;18(4):299-304.

Silveira AM, Lopes HP, Siqueira JF Jr, Macedo SB, Consolaro A.

<u>The clinical efficacy of the local,</u> <u>deep insufflation of an oxygen-</u> <u>ozone mixture in the prevention and</u> <u>treatment of infections in the</u> <u>locomotor system.</u>

Ortop Traumatol Rehabil. 2001;3(4):552-6.

Białoszewski D, Kowalewski M.

Superficially, longer, intermittent ozone theraphy in the treatment of the chronic, infected wounds.

Ortop Traumatol Rehabil. 2003 Oct 30;5(5):652-8.

Białoszewski D, Kowalewski M. <u>Therapeutic efficacy of ozone in</u> <u>patients with diabetic foot.</u>

Eur J Pharmacol. 2005 Oct 31;523(1-3):151-61

Martínez-Sánchez G, Al-Dalain SM, Menéndez S, Re L, Giuliani A, Candelario-Jalil E, Alvarez H, Fernández-Montequín JI, León OS.

Modulation of cutaneous wound healing by ozone: differences between young and aged mice.

Toxicol Lett. 2006 Jan 5;160(2):127-34. Epub 2005 Aug 29.

Lim Y, Phung AD, Corbacho AM, Aung HH, Maioli E, Reznick AZ, Cross CE, Davis PA, Valacchi G. <u>Intravesical ozone therapy for</u> progressive radiation-induced <u>hematuria.</u>

J Altern Complement Med. 2005 Jun;11(3):539-41.

Clavo B, Gutiérrez D, Martín D, Suárez G, Hernández MA, Robaina F. Can the combination of localized "proliferative therapy" with "minor ozonated autohemotherapy" restore the natural healing process?

Med Hypotheses. 2005;65(4):752-9.

Gracer RI, Bocci V.

<u>Major ozonated</u> autohemotherapy in chronic limb ischemia with ulcerations.

J Altern Complement Med. 2005 Apr;11(2):363-7.

de Monte A, van der Zee H, Bocci V. Effects of sodium hypochlorite and ozone on healing of intestinal anastomosis in simulated strangulation colorectal obstruction. Bull Exp Biol Med. 2004 Jan;137(1):103-5.

Lelyanov AD, Sergienko VI, Ivliev NV, Emel'yanov VV, Guseva ED.

Ozone treatment for radiotherapy skin reactions: is there an evidence base for practice?

> Eur J Oncol Nurs. 2002 Dec;6(4):220-7.

Jordan L, Beaver K, Foy S.

[Effect of ozone on antibiotic sensitivity of microorganisms]

Stomatologiia (Mosk). 2003;82(2):36-8

Daulbaeva AA, Baĭzakova GT. Studies on the biological effects of ozone: 10. Release of factors from ozonated human platelets.

Mediators Inflamm. 1999;8(4-5):205-9.

Valacchi G, Bocci V.

[Wound treatment using the flow of an ozonized solution under high pressure]

> Khirurgiia (Mosk). 1998;(8):23-4

Bulynin VI, Ermakova AI, Glukhov AA, Mozhurov IP.

Effects of ozone on how well splitthickness skin grafts according to Thiersch take in war wounds. Results of prospective study.

> Langenbecks Arch Chir. 1995;380(3):144-8.

Turcić J Hancević J, Antoljak T, Zic R, Alfirević I. A physicochemical investigation on the effects of ozone on blood

<u>Travagli, V.a</u>, <u>Zanardi, I.a</u>, <u>Silvietti, A.<sup>b</sup>, Bocci, V.</u> Bocci, V.A.

Scientific and Medical Aspects of Ozone Therapy. State of the Art

(2006) *Archives of Medical Research*, 37 (4), pp. 425-435

### Bocci, V.

Is it true that ozone is always toxic? The end of a dogma

(2006) *Toxicology and Applied Pharmacology*, 216 (3), pp. 493-504. Uppu, R.M., Cueto, R., Squadrito, G.L., Pryor, W.A.

What does ozone react with at the air/lung interface? Model studies using human red blood cell membranes

(1995) Archives of Biochemistry and Biophysics, 319 (1), pp. 257-266. Bocci, V., Luzzi, E., Corradeschi, F., Paulesu, L., Rossi, R., Cardaioli, E., Di Simplicio, P. Studies on the biological effects of

ozone: 4. Cytokine production and glutathione levels in human erythrocytes

(1993) *Journal of Biological Regulators and Homeostatic Agents*, 7 (4), pp. 133-138.

Mudd, J.B., Dawson, P.J., Santrock, J.

Ozone does not react with human erythrocyte membrane lipids

(1997) Archives of Biochemistry and Biophysics, 341 (2), pp. 251-258. Bocci, V., Valacchi, G., Corradeschi, F., Fanetti, G. Studies on the biological effects of ozone: 8. Effects on the total antioxidant status and on interleukin-8 production (1998) *Mediators of Inflammation*, 7 (5), pp. 313-317. Travagli, V., Zanardi, I., Bocci, V. A realistic evaluation of the action of ozone on whole human blood

(2006) International Journal of Biological Macromolecules, 39 (4-5), pp. 317-320

Bocci, V.,

Ozone (2005) *A New Medical Drug*.

Springer, Dordrecht, The Netherlands Biedunkiewicz, B., Lizakowski, S., Tylicki, L., Skiboeska, A., Nieweglowski, T., Chamienia, A., Debska-Slizien, A., (...), Rutkowski, B. Blood Coagulation Unaffected by Ozonated Autohemotherapy in Patients on Maintenance Hemodialysis (2006) Archives of Medical Research, 37 (8), pp. 1034-1037.

Randomised, double-blinded, placebo-controlled, clinical trial of ozone therapy as treatment of sudden sensorineural hearing loss

Bocci, V., Travagli, V., Zanardi, <u>I.</u> Application of medical ozone in endodontic practice

Stomatologiia (Mosk) (Russia 2008;87(6):24-6 (ISSN: 0039-1735)

Bezrukova IV; Petrukhina NB; Dmitrieva NA; Snegirev MV Ozone and its usage in general medicine and dentistry. A review article.

Prague Med Rep 2008;109(1):5-13 (ISSN: 1214-6994) Seidler V; Linetskiy I; Hubalkova H; Stankova H; Smucler R; Mazanek J Effectiveness of ozone against endodontopathogenic microorganisms in a root canal biofilm model. Int Endod J 2009 Jan;42(1):3-13 (ISSN: 1365-2591)

Huth KC; Quirling M; Maier S; Kamereck K; Alkhayer M; Paschos E; Welsch U; Miethke T; Brand K; Hickel R

Reduction by gaseous ozone of Salmonella and microbial flora associated with fresh-cut cantaloupe.

Food Microbiol 2008 Jun;25(4):558-65 (ISSN: 1095-9998)

Selma MV; Ibanez AM; Cantwell M; Suslow T Efficacy of ozonated and electrolyzed oxidative waters to decontaminate hides of cattle before slaughter.

J Food Prot 2005 Jul;68(7):1393-8 (ISSN: 0362-028X) Bosilevac JM; Shackelford SD; Brichta DM; Koohmaraie M Lelyanov AD, Sergienko VI, Ivliev NV, et al. Effects of sodium hypochlorite and ozone on healing of intestinal anastomosis in simulated strangulation colorectal obstruction. Bull Exp Biol Med Jan 2004, 137(1) p103-5

### Kim HS, Noh SU, Han YW, et al.

Therapeutic effects of topical application of **ozone** on acute cutaneous wound **healing**.

J Korean Med Sci, Jun 2009, 24(3) p<u>368-74</u> de Monte A, van der Zee H, Bocci V

Major ozonated autohemotherapy in chronic limb ischemia with ulcerations. J Altern Complement Med (United States), Apr 2005, 11(2) p363-7 Grigor'ian AS, Grigor'iants LA, Guchetl' MN

Experimental-morphological study of the anti-inflammatory action of **ozone**-perfluorane complex application Stomatologiia (Mosk) (Russia (federation)), 2008, 87(2) p4-9 Ozone-initiated disinfection kinetics of Escherichia coli in water.

J Environ Sci Health A Tox Hazard Subst Environ Eng 44;1(48-56 S1093-4529

Zuma F; Lin J; Jonnalagadda SB

Effectiveness of ozonated water on Candida albicans, Enterococcus faecalis, and endotoxins in root canals.

Oral Surg Oral Med Oral Pathol Oral Radiol Endod 2008 Mar;105(3):e85-91 (ISSN: 1528-395X) Cardoso MG; de Oliveira LD; Koga-Ito CY; Jorge AO Li LJ, Yang YG, Zhang ZL, et al.

Protective effects of medical ozone combined with traditional Chinese medicine against chemically-induced hepatic injury in dogs. World J Gastroenterol (China), Dec 7 2007, 13(45) p5989-94

### Thabet SS, Thabet HS, Atalla SS

Efficacy of medical **ozone** in attenuation of murine Schistosomiasis mansoni infection morbidity. J Egypt Soc Parasitol (Egypt), Dec 2007, 37(3) p915-44 Fan L, Song J, McRae KB, et al.

Gaseous ozone treatment inactivates Listeria innocua in vitro. J Appl Microbiol (England), Dec 2007, 103(6) p2657-63 Lin YC, Juan HC, Cheng YC

Ozone exposure in the culture medium inhibits enterovirus 71 virus replication and modulates cytokine production in rhabdomyosarcoma cells. Antiviral Res (Netherlands), Dec 2007, 76(3) p241-51

Balcioglu IA, Tarlan E, Kivilcimdan C, et al.

Merits of **ozonation** and catalytic **ozonation** pre-treatment in the algal treatment of pulp and paper mill effluents. J Environ Manage (England), Dec 2007, 85(4) p918-26

### Bialka KL, Demirci A

Decontamination of Escherichia coli O157:H7 and Salmonella enterica on blueberries using ozone and pulsed UV-light. J Food Sci (United States), Nov 2007, 72(9) pM391-6 Nakada N, Shinohara H, Murata A, et al. Removal of selected pharmaceuticals and personal care products (PPCPs) and endocrinedisrupting chemicals (EDCs) during sand filtration and **ozonation** at a municipal sewage treatment plant. Water Res (England), Nov 2007, 41(19) p4373-82 Rodrigues KL, Cardoso CC, Caputo LR, et al.

Cicatrizing and antimicrobial properties of an **ozonised** oil from sunflower seeds. Inflammopharmacology (Netherlands), 2004, 12(3) p261-70 Silveira AM, Lopes HP, Siqueira JF, et al.

Periradicular repair after two-visit endodontic treatment using two different intracanal medications compared to single-visit endodontic treatment. Braz Dent J (Brazil), 2007, 18(4) p299-304 Rae ID

Ozonised oils as disinfectants.

Ambix (England), Mar 2006, 53(1) p3-20

 Professor Bocci has proven a striking cleansing effect with improved oxygenation and enhanced healing of soreness to diabetic ulcers, burns, traumatic and surgical wounds, abscesses and skin reactions after radiotherapy.

Grootveld M, Silwood C, Lynch E.

Biofactors 2006; 27: 5 – 18.

Oxidative Consumption of Biomolecules by Ozone; Clinical Relevence.

- The application of Ozone in dentistry: A systematic review of literature
- Azarpazhooh and Limeback, J Dent 2008; 36: 104 – 116
- Good evidence of Ozone biocompatibility with human oral epithelial cells, gingival fibroblasts and periodontal cells
- Ozone removes micro-organisms from dental unit water lines, the oral cavity and dentures
- Good evidence of the prophylactic application of Ozone in Restorative Dentistry

### Ebensberger U et al

Use of Ozonated water for avulsed teeth PCNA-expression of cementoblasts and fibroblasts on the root surface after extraoral rinsing with ozonized water for decontamination.

Dent Traumatol 2002; 18: 262 – 266.

- Ozonated Water had no negative effect on periodontal cells remaining on the root surface after irrigation for 2 minutes
- Irrigate all avulsed teeth with Ozonated Water for decontamination and to have a positive effect on the cementoblasts.

J Dent Res 86 (5): 451–456, 2007 Huth KC et al. Effect of Aqueous Ozone on	<ul> <li>Scientific and medical aspects of Ozone therapy. State of the art.</li> </ul>	
the NF-kB System	Professor Velio Bocci	
Aqueous Ozone exerts inhibitory effects on the NF-kB system suggesting it has an anti-inflammatory capacity.	• Arch Med Res 2006; 37: 425 – 435.	The state
IkBa proteolysis, cytokine expression and kB-dependent transcription were prevented.	<ul> <li>Beneficial effects of Ozone discussed</li> </ul>	and the second
		A BORTOS

	<ul> <li>The use of Ozone in dentistry and maxillofacial surgery</li> </ul>	[Preventive use of ozone, short waves, and laser therapy alone and in combination in
Use of Ozonated water for oral surgery	• Stubinger et al	<u>early postoperative period</u> after dental implantation]
	• Quintessence Int 2006, 37 (5):353 – 359.	Vopr Kurortol Fizioter Lech Fiz Kult. 2002 Nov-Dec;(6):17-9
		Korzhachkina NB, Radzievskiĭ SA, Olesova VN.

Clavo B, Gutierrez D, Martin D, et al.	Ozone therapy in the treatment of avascular bisphosphonate – related jaw osteonecrosis	<ul> <li>New therapeutic protocol in the treatment of avascular necrosis of the jaws</li> </ul>
Intravesical <b>ozone</b> therapy for progressive radiation-induced hematuria. J Altern Complement Med (United States), Jun 2005, 11(3) p539-41	<ul> <li>Agrillo et al. J Craniofac Surg 2007; 18: 1071 – 1075</li> <li>Ozone increased the complete healing of lesions with the disappearance of symptoms and brings lesion progression down to zero. Ozone's "benefits were</li> </ul>	<ul> <li>Agrillo et al.</li> <li>J Craniofac Surg 2006; 17: 1080 – 1083</li> </ul>

• J Craniofac Surg. 2007 Sep;18(5):1068-70.	Role of Ozone therapy in the treatment of osteonecrosis of the jaws in multiple myeloma	Ozone accelerates wound healing
<ul> <li>Ozone therapy in extractive surgery on patients treated with bisphosphonates.</li> </ul>	• Petrucci et al	• Chronic leg ulcers. Aust Farm Physician 1985; 14: 292 - 298
<u>Agrillo A</u> <u>Sassano P</u> <u>Rinna C</u> <u>Priore P</u> <u>Iannetti G</u>	• Haematologica 2007; 92: 1289 - 1290	• Skin grafts. Langenbecks Arch Chir 1995; 42: 2 – 6.

<ul> <li>Ozone activates cellular metabolism. Masui 1993; 42: 2 – 6.</li> </ul>	<ul> <li>Wounds with delayed healing were positively</li> </ul>	<ul> <li>Ozonated water clearly accelerates the healing of the human oral mucosa.</li> <li>Professor Filippi Dtsch Zahnarztl Z. 2001;</li> </ul>
Ozone raises intracellular ATP.     Toxicology 1001; 70: 105 202	influenced by Ozonized	56: 104 – 108.
Toxicology 1991, 70, 195 – 202.	water	Ozonated water clearly reduced     complications after impacted third molar
Ozone increases cytokines relevant to wound healing, especially	• Sader et al	surgery. Professor Filippi Dtsch Zahnarztl Z. 1999; 54: 619 - 622.
Transforming Growth Factor. Professor Bocci et al, J Biol Regul	Dtsch Z Mund Kiefer	
Homeost Agents 1994; 8: 108 – 112.	GesichtsChir 1996; 20: 60.	<ul> <li>5 ppm Ozonated Water will eliminate all microorganisms in 3 seconds. Appl</li> </ul>
		Environ Microbiol 1982; 43: 603 – 608.

• The dual action of Ozone on the skin		Water disinfection of dental treatment units using Ozone
• Valacchi et al		• Professor A Filippi et al
<ul> <li>British Journal of Dermatology 2005; 153: 1096 – 1100</li> </ul>	Use of Ozone to treat Dental Unit Water Lines	• Dtsch Zahnarzti Z. 1991; 46 (7): 485 – 487.
<ul> <li>Proven beneficial effect after exposure to Ozone or Ozonated oils to chronic wounds</li> </ul>		<ul> <li>Disinfection proven using Ozonated Water</li> </ul>

• Ozone is the most effective disinfectant for dental treatment units: results after 8 years of comparison

- Professor Filippi
- Ozone Sci Eng 1997; 19: 527

Use of Ozonated water for cleaning Dentures and Impressions Arita M, Nagayoshi M, Fukuizumi T, Okinaga T, Masumi S, Morikawa M, Kakinoki Y, <u>Nishihara T.</u> Microbicidal efficacy of ozonated water against *Candida albicans* adhering to acrylic denture plates. Oral Microbiol Immumol 2005; 20: 206 - 210

 Disinfection of removable dentures using Ozone

- Murakami et al
- Dental Materials 1996;15: 220 225
- Testing of a denture cleaning method using Ozone
- Oizumi et al
- J Med Dent Sci 1998; 45: 135 139

Ozone uses with Caries Professor Sebastian Ciancio, USA Biological Therapies in Dentistry 2005 "Studies with over 2000 Patients have shown remineralisation of Healozone treated permanent and deciduous occlusal surfaces as well as root caries". "The HealOzone is a novel new and painless way to treat early caries".

Amna Al Shamsi

- PhD thesis 2007
- Ozone significantly reduces caries incidence around orthodontic brackets (9% versus 28%)

Dahrhardt J, Jaeggi T, Lussi A. Treating open carious lesions in anxious children with ozone. A prospective controlled clinical study. American Journal of Dentistry 2006; 19: 267 – 270.

Huth K, Paschos E, Brand and Hickel R. Effect of ozone on non cavitated fissure carious lesions in permanent molars. A controlled prospective study.

American Journal of Dentistry 2005; 18: 223 – 270. Effect of ozone treatment on different cariogenic microorganisms in vitro.

Swed Dent J 2008;32(3):139-47 (ISSN: 0347-9994) Fagrell TG; Dietz W; Lingstrom P; Steiniger F; Noren JG The Use of Ozone in Dentistry and Medicine. Part 1.

Baysan A and Lynch E.

Primary Dental Care, 12; 2: April 2005, 47 – 52.

The use of Ozone in Dentistry and Medicine.

• Part 2.

- Baysan and Lynch
- Primary Dental Care 2006; 13: 37 – 41

Anti-microbial effects of a novel ozone generating device on micro-organisms associated with root carious lesions

A. Baysan, R. Whiley and E. Lynch Caries Research 2000;34:498-501.

99% microbial killing achieved after ozone treatment

HealOzone treatment for deep caries or as an alternative to stepwise excavation

Studies in London and Isle of Wight Reversal of root caries using Ozone - A 12 month longitudinal study

A. Baysan (London)

**PhD Thesis University of London 2003** 

### Am J Dent 2004; 17: 54 - 60.

- Baysan A and Lynch E.
- 91% Reversal of Caries at 5.5 months after ozone treatment.
- 99% Reduction in Microorganisms at 5.5 months.

•The Root Caries deeper then 2mms did not reverse. Gerodontology 2003; 20: 106 – 114.

Clinical reversal of root caries using ozone, double blind, randomised, controlled 18-month trial.

- Julian Holmes
- 100% Reversal of Non Cavitated Root Caries at 18 months

Antibacterial effect of ozone on cariogenic bacterial species.

J Dent. 2009 Jun;37(6):449-53

Johansson E Claesson R van Dijken JW

The influence of Healozone on microleakage and fissure penetration of different sealing materials.

Coll Antropol. 2009 Mar;33(1):157-62.

Dukić W Dukić OL Milardović S

Ozonated water improves lipopolysaccharide-induced responses of an odontoblast-like cell line.

J Endod. 2009 May;35(5):668-72

<u>Noguchi F</u><u>Kitamura C</u><u>Nagayoshi M</u> <u>Chen KK</u><u>Terashita M</u><u>Nishihara T</u>

The inability of Streptococcus mutans and Lactobacillus acidophilus to form a biofilm in vitro on dentine pretreated with ozone.

Aust Dent J. 2008 Dec;53(4):349-53

Knight GM McIntyre JM Craig GG Mulyani Zilm PS Treating sensitive cervical areas with ozone. A prospective controlled clinical trial.

Am J Dent. 2008 Apr;21(2):74-6

<u>Dähnhardt JE Gygax M Martignoni B</u> <u>Suter P</u> Lussi A Skaug N, Strand G, Nielsen O

**ORCA & Caries Research 2006** 

85% Reduction in Mutans Streptococci

78.6% Reduction in Lactobacilli

### Khairul Matin, Junji Tagami

99% Microbial killing with HealOzone

Dept of Restorative Sciences and Cariology Tokyo Medical and Dental University Professor Michael Noack and Suzanne Kneist

**ORCA & Caries Research 2006** 

99% reduction in microorganisms following 20 seconds of HealOzone treatment HealOzone treatment for deep caries or as an alternative to stepwise excavation

Studies in London and Isle of Wight

- Dr Dan Mc Kenna
- Proven successful management of Periimplantitis using Ozone 2008 and 2009
- IADR first prize 2009

In vitro reduction of mutans streptococci by means of ozone gas application.

Quintessence Int 2008 Nov;39(10):827-31

Castillo A; Galindo-Moreno P; Avila G; Valderrama M; Liebana J; Baca P Antimicrobial potential of ozone in an ultrasonic cleaning system against Staphylococcus aureus.

Braz Dent J 2006;17(2):134-8

Estrela C; Estrela CR; Decurcio Dde A; Silva JA; Bammann LL

### GROOTVELD M, SILWOOD C AND LYNCH E.

HIGH RESOLUTION NMR INVESTIGATIONS OF THE OXIDATIVE CONSUMPTION OF BIOMOLECULES USING OZONE: RELEVENCE TO THE THERAPEUTIC APPLICATIONS IN CLINICAL DENTISTRY.

**Biofactors 27; 5 - 18, 2006** 

<sup>1</sup>H NMR investigations of the molecular nature of lowmolecular-mass calcium ions in biofluids.

Silwood CL, Grootveld M, Lynch E.

J Biol Inorg Chem. 2002 Jan;7(1-2):46-57 Silwood C, Lynch E, Seddon S, Sheerin A, Claxson A. and Grootveld M. (1999)

<sup>1</sup>H NMR Analysis of Microbial-Derived Organic Acids in Carious Lesions.

NMR in Biomedicine 12: 345-356.



The Journal of Adhesive Dentistry February 2005,	<ul> <li>The impact of Ozone treatment on enamel physical properties</li> </ul>	<ul> <li>The effects of Ozone gas application on Shear bond</li> </ul>
29-32	• Celiberti et al	strength of orthodontic brackets to enamel
Effect of Ozone on Enamel and Dentin Bond Strength	• Am J Dent 2006; 19: 67 – 72.	• Al Shamsi, Cunningham, Lamey and Lynch
PR Schmidlin, Jörg Zimmermann and Andreas Bindl	<ul> <li>No effect of Ozone on enamel physical properties and its effects on sealing ability</li> </ul>	
		• Am J Dent 2008; 21: 35 - 38

### Amna Al Shamsi

- PhD thesis 2007
- Ozone significantly reduces caries incidence around orthodontic brackets (9% versus 28%)

Effect of ozone gas application on the mechanical properties of dental adhesives bonded to dentin.

> Dent Mater. 2008 Oct;24(10):1428-34.

Magni E, Ferrari M, Hickel R, Huth KC, Ilie N. Antibacterial effect of an ozone device and its comparison with two dentin-bonding systems.

Eur J Oral Sci 2006 Aug;114(4):349-53 (ISSN: 0909-8836)

Polydorou O; Pelz K; Hahn P

### Onisor I, Bouillaguet S, Krejci I

Influence of different surface treatments on marginal adaptation in enamel and dentin.

J Adhes Dent (England), Jun 2007, 9(3) p297-303



## Ozone uses with Superbugs

### MRSA and C Difficile

Prions

Ozone uses to clean toothbrushes A quantitative approach to the effectiveness of ozone against microbiota organisms colonizing toothbrushes.

J Dent. 2008 Aug;36(8):600-5. Epub 2008 May 27

<u>Bezirtzoglou E</u> <u>Cretoiu SM</u> <u>Moldoveanu M</u> <u>Alexopoulos A</u> <u>Lazar V</u> <u>Nakou M</u>

Bleaching; Ozone and hydrogen peroxide combinations Effect of ozone, chlorine and hydrogen peroxide on the elimination of colour in treated textile wastewater by MBR.

> Water Sci Technol 2004;49(4):299-303

Brik M; Chamam B; Schoberl P; Braun R; Fuchs W Effect of oxygen, ozone and hydrogen peroxide bleaching stages on the contents and composition of extractives of Eucalyptus globulus kraft pulps. Bioresour Technol 2006 Feb;97(3):420-8

Freire CS; Silvestre AJ; Pascoal Neto C; Evtuguin DV

Dechlorination of chlorophenols found in pulp bleach plant E-1 effluents by advanced oxidation processes.

> Bioresour Technol 2005 May;96(8):897-906

Wang R; Chen CL; Gratzl JS

Rate of dibutylsulfide decomposition by ozonation and the O3/H2O2 advanced oxidation process.

J Hazard Mater 2009 May 30;164(2-3):1364-71

Popiel S; Nalepa T; Dzierzak D; Stankiewicz R; Witkiewicz Z Dehalogenation, degradation and mineralization of diuron by peroxone (peroxide/ozone) treatment.

J Environ Sci Health A Tox Hazard Subst Environ Eng 44;6(630-8 S1093-4529

Catalkaya EC; Kargi F

Kinetics study on photochemical oxidation of polyacrylamide by ozone combined with hydrogen peroxide and ultraviolet radiation. J Environ Sci (China)

2006;18(4):660-4

Ren GM; Sun DZ; Chung JS

Elimination of Listeria monocytogenes biofilms by ozone, chlorine, and hydrogen peroxide.

J Food Prot 2005 Mar;68(3):494-8

Robbins JB; Fisher CW; Moltz AG; Martin SE An ozone/hydrogen peroxide/microwave-enhanced advanced oxidation process for sewage sludge treatment.

J Environ Sci Health A Tox Hazard Subst Environ Eng 42;8(1177-81 S1093-4529

Yin **G**; Liao PH; Lo KV

Oxidative degradation of Nnitrosodimethylamine by conventional ozonation and the advanced oxidation process ozone/hydrogen peroxide.

Water Res 2007 Feb;41(3):581-90

Lee C; Yoon J; Von Gunten U

Treatment of volatile organic chemicals on the EPA Contaminant Candidate List using ozonation and the O3/H2O2 advanced oxidation process. Environ Sci Technol 2006 Apr 15;40(8):2734-9

Chen WR; Sharpless CM; Linden KG; Suffet IH Ozonation with ultrasonic enhancement of p-nitrophenol wastewater.

J Zhejiang Univ Sci B 2005 May;6(5):319-23

Xu XW; Shi HX; Wang DH

Applications of advanced oxidation processes: present and future.

> Water Sci Technol 2004;49(4):227-33

Suty H; De Traversay C; Cost M Degradation of a commercial textile biocide with advanced oxidation processes and ozone.

J Environ Manage 2007 Jan;82(2):145-54 (ISSN: 0301-4797)

**Arslan-Alaton I** 

A comparison of single oxidants versus advanced oxidation processes as chlorine-alternatives for wild blueberry processing (Vaccinium angustifolium). Int J Food Microbiol 2007 May 1;116(1):25-31

Crowe KM; Bushway AA; Bushway RJ; Davis-Dentici K; Hazen RA Seo S, King JM, Prinyawiwatkul W

Simultaneous depolymerization and decolorization of chitosan by **ozone** treatment. J Food Sci (United States), Nov 2007, 72(9) pC522-6 Effects of ozone, ultraviolet and peracetic acid disinfection of a primary-treated municipal effluent on the immune system of rainbow trout (Oncorhynchus mykiss). Comp Biochem Physiol C Toxicol Pharmacol 2008 Aug;148(2):122-7 Hebert N; Gagne F; Cejka P; Bouchard B; Hausler R; Cyr DG; Blaise C; Fournier M

Ozone as Janus: this controversial gas can be either toxic or medically useful.

Mediators Inflamm 2004 Feb;13(1):3-11 (ISSN: 0962-9351) Bocci V

# Ozone gas is an effective and practical antibacterial agent.

Am J Infect Control 2008 Oct;36(8):559-63 (ISSN: 1527-3296) Sharma M; Hudson JB Increase in the ozone decay time in acidic ozone water and its effects on sterilization of biological warfare agents.

J Hazard Mater 2009 Sep 15;168(2-3):1595-601 (ISSN: 1873-3336) Uhm HS; Hong YF; Lee HY; Park YH Application of gaseous ozone for inactivation of Bacillus subtilis spores.

J Air Waste Manag Assoc 2006 Feb;56(2):179-85 (ISSN: 1096-2247) Aydogan A; Gurol MD

Therapeutic effects of topical application of ozone on acute cutaneous wound healing.

J Korean Med Sci. 2009 Jun;24(3):368-74

Kim HS, Noh SU, Han YW, Kim KM, Kang H, Kim HO, Park YM.

[Experimental-morphological study of the anti-inflammatory action of ozone-perfluorane complex application]

Stomatologiia (Mosk). 2008;87(2):4-9 Grigor'ian AS, Grigor'iants LA, Guchetl' MN. The case for oxygenozonetherapy.

Bocci V.

Br J Biomed Sci. 2007;64(1):44-9

Periradicular repair after two-visit endodontic treatment using two different intracanal medications compared to single-visit endodontic treatment.

Braz Dent J. 2007;18(4):299-304.

Silveira AM, Lopes HP, Siqueira JF Jr, Macedo SB, Consolaro A. The clinical efficacy of the local, deep insufflation of an oxygenozone mixture in the prevention and treatment of infections in the locomotor system.

Ortop Traumatol Rehabil. 2001;3(4):552-6.

Białoszewski D, Kowalewski M.

Superficially, longer, intermittent ozone theraphy in the treatment of the chronic, infected wounds.

Ortop Traumatol Rehabil. 2003 Oct 30;5(5):652-8.

Białoszewski D, Kowalewski M.

# Therapeutic efficacy of ozone in patients with diabetic foot.

Eur J Pharmacol. 2005 Oct 31;523(1-3):151-61

Martínez-Sánchez G, Al-Dalain SM, Menéndez S, Re L, Giuliani A, Candelario-Jalil E, Alvarez H, Fernández-Montequín JI, León OS. Modulation of cutaneous wound healing by ozone: differences between young and aged mice.

Toxicol Lett. 2006 Jan 5;160(2):127-34. Epub 2005 Aug 29.

Lim Y, Phung AD, Corbacho AM, Aung HH, Maioli E, Reznick AZ, Cross CE, Davis PA, Valacchi G. <u>Intravesical ozone therapy for</u> progressive radiation-induced <u>hematuria.</u>

J Altern Complement Med. 2005 Jun;11(3):539-41.

Clavo B, Gutiérrez D, Martín D, Suárez G, Hernández MA, Robaina F.

Can the combination of localized "proliferative therapy" with "minor ozonated autohemotherapy" restore the natural healing process?

Med Hypotheses. 2005;65(4):752-9.

Gracer RI, Bocci V.

<u>Major ozonated</u> autohemotherapy in chronic limb ischemia with ulcerations.

J Altern Complement Med. 2005 Apr;11(2):363-7.

de Monte A, van der Zee H, Bocci V. Effects of sodium hypochlorite and ozone on healing of intestinal anastomosis in simulated strangulation colorectal obstruction.

Bull Exp Biol Med. 2004 Jan;137(1):103-5. Lelyanov AD, Sergienko VI, Ivliev NV, Emel'yanov VV, Guseva ED. Ozone treatment for radiotherapy skin reactions: is there an evidence base for practice?

> Eur J Oncol Nurs. 2002 Dec;6(4):220-7.

Jordan L, Beaver K, Foy S.

[Effect of ozone on antibiotic sensitivity of microorganisms]

Stomatologiia (Mosk). 2003;82(2):36-8

Daulbaeva AA, Baĭzakova GT. [Preventive use of ozone, short waves, and laser therapy alone and in combination in early postoperative period after dental implantation]

Vopr Kurortol Fizioter Lech Fiz Kult. 2002 Nov-Dec;(6):17-9

Korzhachkina NB, Radzievskiĭ SA, Olesova VN.

Studies on the biological effects of ozone: 10. Release of factors from ozonated human platelets.

Mediators Inflamm. 1999;8(4-5):205-9.

Valacchi G, Bocci V.

[Wound treatment using the flow of an ozonized solution under high pressure]

> Khirurgiia (Mosk). 1998;(8):23-4

Bulynin VI, Ermakova AI, Glukhov AA, Mozhurov IP. Effects of ozone on how well splitthickness skin grafts according to Thiersch take in war wounds. Results of prospective study.

> Langenbecks Arch Chir. 1995;380(3):144-8.

<u>Turcić J, Hancević J, Antoljak T, Zic</u> <u>R, Alfirević I</u>.

Effect of ozone gas application on the mechanical properties of dental adhesives bonded to dentin.

> Dent Mater. 2008 Oct;24(10):1428-34.

<u>Magni E, Ferrari M, Hickel R,</u> <u>Huth KC, Ilie N</u>. A physicochemical investigation on the effects of ozone on blood

<u>Travagli, V.a</u>, <u>Zanardi, I.a</u>, <u>Silvietti, A.<sup>b</sup>, Bocci, V.</u> Bocci, V.A.

Scientific and Medical Aspects of Ozone Therapy. State of the Art

(2006) *Archives of Medical Research*, 37 (4), pp. 425-435.

### Bocci, V.

Is it true that ozone is always toxic? The end of a dogma

(2006) *Toxicology and Applied Pharmacology*, 216 (3), pp. 493-504. Uppu, R.M., Cueto, R., Squadrito, G.L., Pryor, W.A. What does ozone react with at the air/lung interface? Model studies using human red blood cell membranes

(1995) *Archives of Biochemistry and Biophysics*, 319 (1), pp. 257-266.

Bocci, V., Luzzi, E., Corradeschi, F., Paulesu, L., Rossi, R., Cardaioli, E., Di Simplicio, P. Studies on the biological effects of ozone: 4. Cytokine production and glutathione levels in human erythrocytes (1993) *Journal of Biological Regulators and Homeostatic Agents*, 7 (4), pp. 133-138.

Mudd, J.B., Dawson, P.J., Santrock, J. Ozone does not react with human erythrocyte membrane lipids (1997) Archives of Biochemistry and Biophysics, 341 (2), pp. 251-258.

Bocci, V., Valacchi, G., Corradeschi, F., Fanetti, G. Studies on the biological effects of ozone: 8. Effects on the total antioxidant status and on interleukin-8 production (1998) *Mediators of Inflammation*, 7 (5), pp. 313-317. Travagli, V., Zanardi, I., Bocci, V. A realistic evaluation of the action of ozone on whole human blood (2006) International Journal of Biological Macromolecules, 39 (4-5), pp. 317-320

Bocci, V.,

Ozone (2005) *A New Medical Drug*.

Springer, Dordrecht, The Netherlands Biedunkiewicz, B., Lizakowski, S., Tylicki, L., Skiboeska, A., Nieweglowski, T., Chamienia, A., Debska-Slizien, A., (...), Rutkowski, B. Blood Coagulation Unaffected by Ozonated Autohemotherapy in Patients on Maintenance Hemodialysis (2006) Archives of Medical Research, 37 (8), pp. 1034-1037.

Randomised, doubleblinded, placebocontrolled, clinical trial of ozone therapy as treatment of sudden sensorineural hearing loss <u>Bocci, V., Travagli, V.,</u> <u>Zanardi, I.</u> Dehalogenation, degradation and mineralization of diuron by peroxone (peroxide/ozone) treatment.

J Environ Sci Health A Tox Hazard Subst Environ Eng 44;6(630-8 S1093-4529

Catalkaya EC; Kargi F

Effect of oxygen, ozone and hydrogen peroxide bleaching stages on the contents and composition of extractives of Eucalyptus globulus kraft pulps. Bioresour Technol 2006 Feb;97(3):420-8

Freire CS; Silvestre AJ; Pascoal Neto C; Evtuguin DV Kinetics study on photochemical oxidation of polyacrylamide by ozone combined with hydrogen peroxide and ultraviolet radiation. J Environ Sci (China)

2006;18(4):660-4

Ren GM; Sun DZ; Chung JS

Elimination of Listeria monocytogenes biofilms by ozone, chlorine, and hydrogen peroxide.

J Food Prot 2005 Mar;68(3):494-8

Robbins JB; Fisher CW; Moltz AG; Martin SE An ozone/hydrogen peroxide/microwave-enhanced advanced oxidation process for sewage sludge treatment.

J Environ Sci Health A Tox Hazard Subst Environ Eng 42;8(1177-81 S1093-4529

Yin G; Liao PH; Lo KV

Oxidative degradation of Nnitrosodimethylamine by conventional ozonation and the advanced oxidation process ozone/hydrogen peroxide.

Water Res 2007 Feb;41(3):581-90

Lee C; Yoon J; Von Gunten U

Treatment of volatile organic chemicals on the EPA Contaminant Candidate List using ozonation and the O3/H2O2 advanced oxidation process. Environ Sci Technol 2006 Apr 15;40(8):2734-9

Chen WR; Sharpless CM; Linden KG; Suffet IH Dechlorination of chlorophenols found in pulp bleach plant E-1 effluents by advanced oxidation processes. Bioresour Technol 2005 May;96(8):897-906

Wang R; Chen CL; Gratzl JS

Rate of dibutylsulfide decomposition by ozonation and the O3/H2O2 advanced oxidation process.

J Hazard Mater 2009 May 30;164(2-3):1364-71

Popiel S; Nalepa T; Dzierzak D; Stankiewicz R; Witkiewicz Z Ozonation with ultrasonic enhancement of p-nitrophenol wastewater.

J Zhejiang Univ Sci B 2005 May;6(5):319-23

Xu XW; Shi HX; Wang DH

Effect of ozone, chlorine and hydrogen peroxide on the elimination of colour in treated textile wastewater by MBR.

> Water Sci Technol 2004;49(4):299-303

Brik M; Chamam B; Schoberl P; Braun R; Fuchs W Applications of advanced oxidation processes: present and future.

> Water Sci Technol 2004;49(4):227-33

Suty H; De Traversay C; Cost M

In vitro reduction of mutans streptococci by means of ozone gas application.

Quintessence Int 2008 Nov;39(10):827-31

Castillo A; Galindo-Moreno P; Avila G; Valderrama M; Liebana J; Baca P Antimicrobial potential of ozone in an ultrasonic cleaning system against Staphylococcus aureus.

Braz Dent J 2006;17(2):134-8

Estrela C; Estrela CR; Decurcio Dde A; Silva JA; Bammann LL Application of medical ozone in endodontic practice Stomatologiia (Mosk) (Russia 2008;87(6):24-6 (ISSN: 0039-1735) Bezrukova IV; Petrukhina NB; Dmitrieva NA; Snegirev MV

Ozone and its usage in general medicine and dentistry. A review article. Prague Med Rep 2008;109(1):5-13 (ISSN: 1214-6994) Seidler V; Linetskiy I; Hubalkova H; Stankova H; Smucler R; Mazanek J

Antibacterial effect of an ozone device and its comparison with two dentin-bonding systems. Eur J Oral Sci 2006 Aug;114(4):349-53 (ISSN: 0909-8836) Polydorou O; Pelz K; Hahn P Effectiveness of ozone against endodontopathogenic microorganisms in a root canal biofilm model. Int Endod J 2009 Jan;42(1):3-13 (ISSN: 1365-2591) Huth KC; Quirling M; Maier S; Kamereck K; Alkhayer M; Paschos E; Welsch U; Miethke T; Brand K; Hickel R Reduction by gaseous ozone of Salmonella and microbial flora associated with fresh-cut cantaloupe. Food Microbiol 2008 Jun;25(4):558-65 (ISSN: 1095-9998) Selma MV; Ibanez AM; Cantwell M; Suslow T

Ozone therapy in medicine and dentistry. J Contemp Dent Pract 2008;9(4):75-84 (ISSN: 1526-3711) Nogales CG; Ferrari PH; Kantorovich EO; Lage-Marques JL Degradation of a commercial textile biocide with advanced oxidation processes and ozone. J Environ Manage 2007 Jan;82(2):145-54 (ISSN: 0301-4797) Arslan-Alaton I

Ozone-initiated disinfection kinetics of Escherichia coli in water. J Environ Sci Health A Tox Hazard Subst Environ Eng 44;1(48-56 S1093-4529 Zuma F; Lin J; Jonnalagadda SB Effectiveness of ozonated water on Candida albicans, Enterococcus faecalis, and endotoxins in root canals. Oral Surg Oral Med Oral Pathol Oral Radiol Endod 2008 Mar;105(3):e85-91 (ISSN: 1528-395X) Cardoso MG; de Oliveira LD; Koga-Ito CY; Jorge AO

A comparison of single oxidants versus advanced oxidation processes as chlorine-alternatives for wild blueberry processing (Vaccinium angustifolium). Int J Food Microbiol 2007 May 1;116(1):25-31 (ISSN: 0168-1605) Crowe KM; Bushway AA; Bushway RJ; Davis-Dentici K; Hazen RA

Efficacy of ozonated and electrolyzed oxidative waters to decontaminate hides of cattle before slaughter. J Food Prot 2005 Jul;68(7):1393-8 (ISSN: 0362-028X) Bosilevac JM; Shackelford SD; Brichta DM; Koohmaraie M Lelyanov AD, Sergienko VI, Ivliev NV, et al. Effects of sodium hypochlorite and ozone on healing of intestinal anastomosis in simulated strangulation colorectal obstruction. Bull Exp Biol Med Jan 2004, 137(1) p103-5

Kim HS, Noh SU, Han YW, et al. Therapeutic effects of topical application of **ozone** on acute cutaneous wound **healing**. J Korean Med Sci, Jun 2009, 24(3) p368-74 de Monte A, van der Zee H, Bocci V Major ozonated autohemotherapy in chronic limb ischemia with ulcerations. J Altern Complement Med (United States), Apr 2005, 11(2) p363-7 Grigor'ian AS, Grigor'iants LA, Guchetl' MN [Experimentalmorphological study of the anti-inflammatory action of ozone-perfluorane complex application] Stomatologiia (Mosk) (Russia (federation)), 2008, 87(2) p4-9

Clavo B, Gutierrez D, Martin D, et al. Intravesical **ozone** therapy for progressive radiationinduced hematuria. J Altern Complement Med (United States), Jun 2005, 11(3) p539-41

Onisor I, Bouillaguet S, Krejci I Influence of different surface treatments on marginal adaptation in enamel and dentin. J Adhes Dent (England), Jun 2007, 9(3) p297-303

Rodrigues KL, Cardoso CC, Caputo LR, et al. Cicatrizing and antimicrobial properties of an **ozonised** oil from sunflower seeds. Inflammopharmacology (Netherlands), 2004, 12(3) p261-70 Silveira AM, Lopes HP, Siqueira JF, et al. Periradicular repair after twovisit endodontic treatment using two different intracanal medications compared to single-visit endodontic treatment. Braz Dent J (Brazil), 2007, 18(4) p299-304

Rae ID Ozonised oils as disinfectants. Ambix (England), Mar 2006, 53(1) p3-20 Li LJ, Yang YG, Zhang ZL, et al. Protective effects of medical ozone combined with traditional Chinese medicine against chemically-induced hepatic injury in dogs. World J Gastroenterol (China), Dec 7 2007, 13(45) p5989-94 Thabet SS, Thabet HS, Atalla SS Efficacy of medical **ozone** in attenuation of murine Schistosomiasis mansoni infection morbidity. J Egypt Soc Parasitol (Egypt), Dec 2007, 37(3) p915-44 Fan L, Song J, McRae KB, et al. Gaseous **ozone** treatment inactivates Listeria innocua in vitro. J Appl Microbiol (England), Dec 2007, 103(6) p2657-63 Lin YC, Juan HC, Cheng YC Ozone exposure in the culture medium inhibits enterovirus 71 virus replication and modulates cytokine production in rhabdomyosarcoma cells. Antiviral Res (Netherlands), Dec 2007, 76(3) p241-51

Balcioglu IA, Tarlan E, Kivilcimdan C, et al. Merits of ozonation and catalytic ozonation pretreatment in the algal treatment of pulp and paper mill effluents. J Environ Manage (England), Dec 2007, 85(4) p918-26

Bialka KL, Demirci A Decontamination of Escherichia coli O157:H7 and Salmonella enterica on blueberries using **ozone** and pulsed UV-light. J Food Sci (United States), Nov 2007, 72(9) pM391-6 Seo S, King JM, Prinyawiwatkul W Simultaneous depolymerization and decolorization of chitosan by ozone treatment. J Food Sci (United States), Nov 2007, 72(9) pC522-6 Nakada N, Shinohara H, Murata A, et al. Removal of selected pharmaceuticals and personal care products (PPCPs) and endocrine-disrupting chemicals (EDCs) during sand filtration and **ozonation** at a municipal sewage treatment plant. Water Res (England), Nov 2007, 41(19) p4373-82

Effect of ozone treatment on different cariogenic microorganisms in vitro. Swed Dent J 2008;32(3):139-47 (ISSN: 0347-9994) Fagrell TG; Dietz W; Lingstrom P; Steiniger F; Noren JG Ozone as Janus: this controversial gas can be either toxic or medically useful. Mediators Inflamm 2004 Feb;13(1):3-11 (ISSN: 0962-9351) Bocci V

Ozone gas is an effective and practical antibacterial agent. Am J Infect Control 2008 Oct;36(8):559-63 (ISSN: 1527-3296) Sharma M; Hudson JB Effects of ozone, ultraviolet and peracetic acid disinfection of a primarytreated municipal effluent on the immune system of rainbow trout (Oncorhynchus mykiss). Comp Biochem Physiol C Toxicol Pharmacol 2008 Aug;148(2):122-7 (ISSN: 1532-0456) Hebert N; Gagne F; Cejka P; Bouchard B; Hausler R; Cyr DG; Blaise C; Fournier M

Increase in the ozone decay time in acidic ozone water and its effects on sterilization of biological warfare agents. J Hazard Mater 2009 Sep 15;168(2-3):1595-601 (ISSN: 1873-3336) Uhm HS; Hong YF; Lee HY; Park YH

Application of gaseous ozone for inactivation of Bacillus subtilis spores. J Air Waste Manag Assoc 2006 Feb;56(2):179-85 (ISSN: 1096-2247) Aydogan A; Gurol MD

Efficacy of calcium hydroxide, Er:YAG laser or gaseous ozone against Enterococcus faecalis in root canals.

Am J Dent. 2009 Feb;22(1):14-8. Noetzel J Nonhoff J Bitter K Wagner J Neumann K Kielbassa AM Antibacterial effect of ozone on cariogenic bacterial species.

J Dent. 2009 Jun;37(6):449-53

Johansson E Claesson R van Dijken JW The influence of Healozone on microleakage and fissure penetration of different sealing materials.

> Coll Antropol. 2009 Mar;33(1):157-62.

<u>Dukić W</u> <u>Dukić OL</u> <u>Milardović</u> <u>S</u>

Ozonated water improves lipopolysaccharide-induced responses of an odontoblastlike cell line.

J Endod. 2009 May;35(5):668-72

Noguchi F Kitamura C Nagayoshi M Chen KK Terashita M Nishihara T Bactericidal effect of KTP laser irradiation against Enterococcus faecalis compared with gaseous ozone: an ex vivo study.

Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 May;107(5):e73-9

Kuştarci A Sümer Z Altunbaş D Koşum S The inability of Streptococcus mutans and Lactobacillus acidophilus to form a biofilm in vitro on dentine pretreated with ozone.

> Aust Dent J. 2008 Dec;53(4):349-53

Knight GM McIntyre JM Craig GG Mulyani Zilm PS Treating sensitive cervical areas with ozone. A prospective controlled clinical trial.

Am J Dent. 2008 Apr;21(2):74-6

Dähnhardt JE Gygax M Martignoni B Suter P Lussi A

# Therapeutic effects of topical application of ozone on acute cutaneous wound healing.

J Korean Med Sci. 2009 Jun;24(3):368-74. Epub 2009 Jun 12

Kim HS Noh SU Han YW Kim KM Kang H Kim HO Park YM Superficially, longer, intermittent ozone therapy in the treatment of the chronic, infected wounds.

Ortop Traumatol Rehabil. 2003 Oct 30;5(5):652-8.

Białoszewski D Kowalewski M

[New therapeutic strategies for the treatment of difficult wounds]

G Chir. 2008 May;29(5):212-20

<u>Onesti MG</u> <u>Bitonti A</u> <u>Fino P</u> <u>Ciotti M</u> <u>Scuderi N</u> A quantitative approach to the effectiveness of ozone against microbiota organisms colonizing toothbrushes.

J Dent. 2008 Aug;36(8):600-5. Epub 2008 May 27

Bezirtzoglou E Cretoiu SM Moldoveanu M Alexopoulos A Lazar V Nakou M Efficacy of calcium hydroxide, Er:YAG laser or gaseous ozone against Enterococcus faecalis in root canals.

Am J Dent. 2009 Feb;22(1):14-8.

Beneficial effects of pro- /antioxidant-based	Featherstone JD.	• Brostek A
nutraceuticals in the skin rejuvenation techniques.	<ul> <li>The caries balance: the basis for caries management by risk assessment. Oral Health Prev</li> </ul>	• • Early diagnosis and minimally invasive treatment of occlusal
Cell Mol Biol (Noisy-le-grand). 2007 Apr 15;53(1):94-101.	Dent. 2004;2 Suppl 1:259-64.	carles-a clinical approach. Oral Health Prev Dent. 2004;2 Suppl 1:313-9.
<u>de Luca C, Deeva I,</u> <u>Mikhal'Chik E, Korkina L</u> .	<ul> <li>a change in how we manage caries is long over due. The days of drilln'fill are numbered</li> </ul>	•accurate diagnosis more difficult.

# Ozone uses to clean toothbrushes

- A quantitative approach to the effectiveness of ozone against microbiota organisms colonizing toothbrushes.
- J Dent. 2008 Aug;36(8):600-5. Epub 2008 May 27

Bezirtzoglou E, Cretoiu SM, Moldoveanu M, Alexopoulos A, Lazar V, Nakou M. The measurement of root caries for research purposes

Edward Lynch

Journal of Dental Research, 65, 510, 207, 1986.

Beighton D, Lynch E and Heath M (1993)	Lynch E. and Beighton D. (1994)	Lynch E. (1996)
A microbiological study of primary root caries with different treatment needs <i>Journal of Dental Research</i> 73: 623-629 This is the only validated severity index	A comparison of primary root caries lesions classified according to colour.	Relationships between clinical criteria and microflora of primary root caries. In: Early Detection of Dental Caries. Ed. Stookey G.K., 195-243
for root caries		

### Lynch E. (1996)

Antimicrobial management of primary root carious lesions.

Gerodontology 13: 118-129

A pharmaceutical approach to the management of root caries

E. Lynch and A. Baysan

Tissue Preservation and Caries Treatment Quintessence Book 2001, Chapter 3, 81-104.

## LYNCH E et al

Effectiveness of two fluoride dentifrices to arrest root carious lesions. Am J Dent. 2000;13:218-220 LYNCH E, BAYSAN A. Reversal of primary root caries using a dentifrice with a high fluoride content. Caries Res 2001: 35;1:60-64 Lynch E, Grootveld M et al. (1999)

Multicomponent evaluations of the oxidising actions and status of a peroxoborate system using high resolution proton NMR spectroscopy.

Journal of Inorganic Biochemistry 73: 65-84.

### Lynch E, Grootveld M et al. (1997)

Multicomponent spectroscopic investigations of antioxidant consumption by an oral preparation containing the stable free radical species chlorine dioxide (CIO<sub>2</sub>•).

Free Radical Research 26:209-237.

Ozone therapy in the treatment of avascular bisphosphonate-related jaw osteonecrosis.

J Craniofac Surg. 2007 Sep;18(5):1071-5

Agrillo A Ungari C Filiaci F Priore P lannetti G

# Role of ozone therapy in the treatment of osteonecrosis of the jaws in multiple myeloma patients.

Haematologica. 2007 Sep;92(

Petrucci MT Gallucci C Agrillo A Mustazza MC Foà R

### Antimicrobial management

Lynch E. et al. (1997)

Multicomponent spectroscopic investigations of salivary antioxidant consumption by an oral rinse preparation containing the stable free radical species chlorine dioxide (ClO<sub>2</sub>\*).

Free Radical Research 26:209-237.

#### Antimicrobial management

Lynch E . *et al.* (1999)

Multicomponent evaluations of the oxidising actions and status of a peroxoborate-containing tooth-whitening system in whole human saliva using high resolution proton NMR spectroscopy.

*Journal of Inorganic Biochemistry* 73: 65-84

#### Antimicrobial management

Silwood C.J.L., Lynch E., Seddon S., Sheerin A., Claxson A.W.D. And Grootveld M. (1999)

<sup>1</sup>H NMR Analysis of Microbial-Derived Organic Acids in Primary Root Carious Lesions and Saliva.

NMR in Biomedicine 12: 345-356.

### <sup>1</sup>H NMR spectra of root caries



Queen's Univer

### A Pharmaceutical Approach to the Management of Root Caries

### **Edward Lynch**

Professor of Restorative Dentistry and Gerodontology School of Clinical Dentistry Queen's University Belfast



The measurement of root caries for research purposes

Lynch E.

J Dent Res 65, 510, 207, 1986.

Beighton D., Lynch E. and Heath M.R. (1993)

A microbiological study of primary root caries lesions with different treatment needs.

Journal of Dental Research 73: 623-629.

This is the only validated severity index for root caries

Lynch E. and Beighton D. (1993)

Relationships between Mutans streptococci and perceived treatment needs of primary root caries lesions.

Gerodontology 10: 98-104.

Lynch E. and Beighton D. (1994)

A comparison of primary root caries lesions classified according to colour.

Caries Research 28: 233-239.

Lynch E. (1996)

Relationships between clinical criteria and microflora of primary root caries.

In: Early Detection of Dental Caries. Ed. Stookey G.K., 195-243.

Lynch E. (1996)

Antimicrobial management of primary root carious lesions.

Gerodontology 13: 118-129.

A pharmaceutical approach to the management of root caries

E. Lynch and A. Baysan Tissue Preservation and Caries Treatment Quintessence Book 2001, Chapter 3, p. 81-104. Management of root caries using a dentifrice with a high FLUORIDE content
Reversal of primary root caries using dentifrices containing 5,000 and 1,100 ppm fluoride.

Baysan A, Lynch E, Ellwood R, Davies R, Petersson L, Borsboom P. Caries Res. 2001;35:41-46.

Reversal of primary root caries using dentifrices containing 5,000 and 1,100 ppm fluoride - 3 month follow-up

Lynch E, Baysan A, Ellwood R, Davies R, Petersson L, Borsboom P. Amer J Dent 2000;13:218-221.

Management of primary root caries using a dentifrice with a high fluoride content

A. Baysan and E. Lynch **Tissue Preservation and Caries Treatment** Quintessenz Book 2001, Chapter 2, p. 37-48.

A pharmaceutical approach to the management of root caries

E. Lynch and A. Baysan **Tissue Preservation and Caries Treatment** Quintessenz Book 2001, Chapter 3, p. 81-104.

Reversal of primary root caries using a dentifrice with a high fluoride content

E. Lynch and A. Baysan Caries Research





Root caries remineralisation using 5,000 and 1,100 ppm fluoride dentifrices

A. BAYSAN\* and E. LYNCH<sup>1</sup>

The Morita Investigator Award for the Best Clinical cience Presentation from Geriatric Oral Research Grou at the International Association of Dental Research in April 2000 Scie

### Conclusion

The use of dentifrices containing either 5,000 or 1,100 ppm fluoride was associated with the reversal of some of PRCLs.

The use of a dentifrice with a high fluoride content was significantly better to reverse leathery lesions than an 1,100 ppm fluoride dentifrice within 6 months. **Treatment of cervical sensitivity** with a root sealant

A. BAYSAN<sup>1,2\*</sup>, E. LYNCH<sup>1</sup>, S. BRAILSFORD<sup>3</sup> and D. BEIGHTON<sup>3</sup> <sup>1</sup>Restorative Dentistry and Gerodontology, Queen's University Belfast, Northern Ireland, <sup>2</sup>Department of Adult Oral Health, St. Bar and the Royal London School of Medicine and Dentistry, Londor UK, <sup>3</sup>Oral Microbiology, Guy's, King's and St. Thomas Dental Institute, London sity Rart's

Dentine sensitivity is one of the most painful and least predictably treated clinical conditions.

It has been established that dentinal hypersensitivity affects 1 in 6 people. Incidence tends to peak around the third decade of life and is equally divided between men and women.

# Aim

The aim of this study was to assess a new protective root sealant for the treatment of cervical sensitivity.

### **Results**

Sensitivity scores at baseline and at time points of 3, 6 and 19 months



## Conclusions

There was a significant reduction in sensitivity scores compared to baseline after 19 months.

The protective sealant was found to be capable of covering the cervical surface to prevent further wear.

In addition, there was a significant reduction of some representative caries associated micro-organisms in the overlying plaque.

### Understanding the Results \_\_\_\_\_ (5/7)

51 - 90

LED: 3 vellows (optional)

AUDIO: 2 beeps

LCD:

OBSERVATION: Probability of significant carious change beneath the enamel surface warranting specific preventive care.

ADVICE: Preventive Care Advised (PCA) Strongly consider fluoride varnish or (on pit & fissure sites) pit & fissure sealants, with localized oral hygiene, diet advice and personalised review with monitoring at shorter intervals.



Is Ozone useful for Periodontal –Endodontic Lesions?

Cleanability of dental instruments--implications of residual protein and risks from Creutzfeldt-Jakob disease.

Walker JT, Dickinson J, Sutton JM, Raven ND, Marsh PD. Br Dent J. 2007 Oct 13;203(7):395-401.

### WWW.UKDENT.COM

LIMITED TIME SPECIAL OFFER ON NEW HEALOZONE x4 Is Ozone useful for Periodontal –Endodontic Lesions?





### The CarieScan PRO



#### How the PRO works



#### A small electrical current is passed through the tooth and the response is measured. Comparison of the response with the applied signal leads to the impedance.

### nce is an electrical parameter that is influenced by the physical es of the material being studied.

#### **User Guide**



### **Treatment options**

- Dietary Advice
- Oral Hygiene Advice
- Re-mineralising agents Great Oral Health products
- Pit & Fissure Sealants
- Topical Fluorides and Varnishes
- Other agents

### RemoteView

- Wireless ca
- Monitoring featu
- Alternative dis
- Included in system price







### E.LYNCH@WARWICK.AC.UK

THE UNIVERSITY OF

Conclusion

HealOzone is safe





