CLEARFIL™ DC CORE PLUS

ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION

CLEARFIL DC CORE PLUS is a dual-cured (light-cured with self-curing property), radiopaque two-component core build-up material supplied in an automix delivery system. CLEARFIL DC CORE PLUS is available in two shades: Dentin and White.

II. INDICATIONS

Post cementation and core build-up

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

IV. INCOMPATIBILITIES

Do not use eugenol-containing materials for temporary sealing since the eugenol can retard the curing process.

V. PRECAUTIONS

1. Safety precautions

- The product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
- If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue the use of the product and seek medical attention.
- Avoid direct contact with the skin and/or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- If the product comes in contact with human body tissues, take the following actions:
 - < If the product gets in the eye > Immediately wash the eye with copious amounts of water and consult a physician.
 - < If the product comes in contact with the skin or the oral mucosa > Immediately wipe the area with a cotton pledget or a gauze moistened with alcohol; rinse with copious amounts of water.
- Use caution when using the product to prevent swallowing.
- Avoid looking directly at the dental curing light when curing the product.
- Do not reuse the mixing tip and the guide tip to prevent cross-contamination. The mixing tip and the guide tip are single use only. Discard them after use.

2. Handling and manipulation precautions

[CLEARFIL DC CORE PLUS]

- Only CLEARFIL TRI-S BOND PLUS or dual-cured bonding agents can be used with the product. Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS. In the areas where the light cannot reach easily (e.g. adherent areas masked by the post), "CLEARFIL TRI-S BOND PLUS" can be self-cured in contact with the product.
- The product should not be used for purposes other than specified in [III. INDICATIONS].
- The use of the product is restricted to licensed dentist.
- Provide moisture and contamination control with the use of a rubber dam.
- Amalgam, temporary sealing material, lining material etc. remaining in the root canal will prevent the passage of light and the polymerization of the product. Completely remove any lining material when preparing the root canal.
- Clean the root canal sufficiently to prevent poor bonding. If the adherent surface is contaminated with saliva or tissue exudates, wash it thoroughly and dry before cementation or bonding.
- The product should be brought to room temperature for 15 minutes or more after it is removed from the refrigerator. If the syringe is not brought up to room temperature, you will find it difficult to squeeze out the paste.
- Do not use a lentulo spiral to insert the paste into the root canal; this can accelerate the polymerization of the paste beyond the desirable limit.
- Do not use the product in conjunction with other composite resins for core build-up. Intermixing materials may cause a change in physical properties, including a possible decrease, from the properties expected.
- Be careful to prevent unnecessary exposure to direct sunlight or surgery operating lights, otherwise the paste inside the tip may harden, leading to a shortened working time.
- Although the product is a dual-cured core build-up material, it must be light-cured according to the instructions.
- When light-curing the product, note the light-curing depth in the Instructions for Use.
- While putting the product into the root canal, be careful to avoid cross-contamination. Cover the entire syringe with a disposable barrier (e.g. a poly bag) to prevent saliva and blood contamination. Disinfect the syringe by wiping them with an absorbent cotton with alcohol both before and after use.
- Do not immerse the syringe in a disinfectant solution.
- Be careful not to cut your fingers on the sharp edges of instruments.

[Dental light-curing unit]

- Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- The product must be refrigerated (2-8°C / 36-46°F) when not in use, and should be brought to room temperature before using.
- Keep away from extreme heat or direct sunlight.
- The product must be stored in proper places where only dental practitioners can access it.

VI. COMPONENTS

Please see the outside of the package for contents and quantity.

- A Paste and B Paste -Dentin or White
 - Principal ingredients
 - (1) A Paste
 - Bisphenol A diglycidylmethacrylate (Bis-GMA)
 - Hydrophobic aliphatic dimethacrylate
 - Hydrophilic aliphatic dimethacrylate
 - Hydrophobic aromatic dimethacrylate
 - Silanated barium glass filler
 - Silanated colloidal silica
 - Colloidal silica
 - di-Camphorquinone
 - Initiators
 - Pigments
 - (2) B Paste
 - Triéthylène glycol diméthacrylate
 - Hydrophilic aliphatic diméthacrylate
 - Hydrophobic aromatic diméthacrylate
 - Silanated barium glass filler
 - Silanated colloidal silica
 - Aluminium oxide filler
 - Accelerators


The total amount of inorganic filler is approx. 52 vol%. The particle size of inorganic fillers ranges from 0.01 µm to 20 µm.

- Accessories
 - Mixing tip (Mixing tip)
 - Guide tip (L)
 - Guide tip (S)

VII. CLINICAL PROCEDURES

SYRINGE PREPARATION


The device components.


- Turn the cap 1/4 of a turn counterclockwise to align the projections on the cap with the grooves (A) on the syringe. Holding the base of the cap, remove the cap by twisting and pressing downward (B).


- Before each use, extrude small amounts of the two pastes, making sure equal amounts are being dispensed through the two outlets of the syringe, and discard them. If equal amounts of paste are not used, there is a possibility of poor polymerization.


- Attach a mixing tip to the syringe, aligning the projections on it with the grooves in the syringe, as closely as possible.


- Turn the mixing tip 1/4 of a turn clockwise to lock in place.

- Insert a guide tip into the mixing tip until it clicks.

- The mixed pastes will be dispensed through the guide tip when the plunger is pressed.

[NOTE]

- After use, the syringe should be stored with the cap. When you put the cap back on the syringe before storage, make sure the cap is free of paste.
- When replacing an old mixing tip with a new one, turn the mixing tip 1/4 of a turn counterclockwise to align the projections of the mixing tip with the grooves on the syringe. Remove the mixing tip from the syringe by twisting and pressing downward.
- If the paste has hardened, making it difficult to squeeze the paste mixture out of the syringe, remove it using an appropriate instrument.
- When attaching the guide tip to the mixing tip or rotating the guide tip after attaching, use caution and hold the root of the guide tip to avoid bending.

POST CEMENTATION AND CORE BUILD-UP

- Moisture control
 - In order to produce optimal results, avoid contamination of the treatment area with saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.
- Preparing root canal
 - Prepare and clean the root canal opening in the usual manner.
- Try-in of the post
 - A post of the proper diameter is fitted to the prepared root canal and the length of the post adjusted.
- Treatment of post surface
 - When using glass fiber pastes
 - Bisphenol A diglycidylmethacrylate (e.g. CLEARFIL CERAMIC PRIMER or the mixture of CLEARFIL PORCELAIN BOND ACTIVATOR and CLEARFIL TRI-S BOND PLUS).
 - When using metal posts
 - Treat the post surface according to the metal adhesive primer's instructions (e.g. ALLOY PRIMER).
- Syringe preparation
 - Select the desired shade of the product and prepare the syringe and accessories according to "SYRINGE PREPARATION". CLEARFIL DISPENSER (10ml, 1:1) can be used for easy dispensing. For detailed descriptions, refer to the Instructions for Use supplied with the CLEARFIL DISPENSER.
- Post cementation
 - 6-a. When the following adhesives are used (Post cementation by the product)
 - Tooth surface treatment and bonding should be performed according to the instructions for use of CLEARFIL TRI-S BOND PLUS or dual-cured bonding agents (e.g. CLEARFIL DC BOND, CLEARFIL LINER BOND 2V or CLEARFIL PHOTO BOND).
 - CAUTION
 - Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS.
 - Carefully squeeze the paste from the syringe directly into the root canal, being sure to avoid entrapment of air. You must start step 3) within 1 minute after application of the paste, since self-cure of the product is accelerated on contact with a surface on which the bonding agent has been applied.

[CAUTION]

Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS.

[CAUTION]

Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS.

[CAUTION]

Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS.

[CAUTION]

Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS.

[CAUTION]

Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS.

[CAUTION]

Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS.

- Timing for CLEARFIL DC CORE PLUS:
 - < Working time (23°C / 73°F) : 3 minutes >
 - < Setting time (37°C / 99°F) : 6 minutes >
- While putting the product into the root canal, be careful to avoid cross-contamination. Cover the entire syringe with a disposable barrier (e.g. a poly bag) to prevent saliva and blood contamination. Disinfect the syringe by wiping them with an absorbent cotton with alcohol both before and after use.

- Insert the post into the canal.

- Cure the paste using the dental curing unit according to the table below. It is advisable to light-cure the margins well for strong bond.


Table: Relationship between curing time and curing depth for each dental curing unit.

Type	Curing time	Curing depth
Halogen	10 seconds	1.5 mm
	20 seconds	2.0 mm
LED	10 seconds	1.5 mm
	20 seconds	2.0 mm

*Dental curing unit

Type	Light source	Wavelength range and light intensity
Halogen	Halogen lamp	Light intensity ²⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm
LED	Blue LED ¹⁾	Light intensity ³⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm

- Peak of emission spectrum: 450 - 480 nm
- Evaluated according to ISO 10650-1.
- Evaluated according to ISO 10650-2.

6-b. When PANAVIA F 2.0, CLEARFIL SA CEMENT or CLEARFIL ESTHETIC CEMENT EX is used.

Cement the post into the root canal according to the instructions for use of PANAVIA F 2.0, CLEARFIL SA CEMENT or CLEARFIL ESTHETIC CEMENT EX.

6-c. When another luting cement is used.

- Cement the post into the root canal according to the instructions for use of the luting cement.
- Apply the bonding agent to the entire adherent surface according to the Instructions for Use of the bonding system (e.g. CLEARFIL TRI-S BOND PLUS, CLEARFIL DC BOND, CLEARFIL LINER BOND 2V or CLEARFIL PHOTO BOND).

7. Core build-up

- After cementing the post in the root canal, squeeze the paste from the syringe directly around the post.


[CAUTION]

While squeezing the paste around the post, be careful to avoid cross-contamination. Cover the entire syringe with a disposable barrier (e.g. a poly bag) to prevent saliva and blood contamination. Disinfect the syringe by wiping them with an absorbent cotton with alcohol both before and after use.

[NOTE]

If necessary, a commercial matrix can be used to stabilize PASTE or tack curing with light for 5 seconds may be used to perform the initial polymerization of incremental paste placement.

- After placing the paste around the post, light-cure the paste from both the lingual and the labial (buccal) sides, using a dental curing unit. See table "Relationship between curing time and curing depth for each dental curing unit" in 6-a. Leave the paste more than 6 minutes after light-curing when the thickness of the paste is greater than the specified curing depth on the table.

- Preparation of an abutment tooth
 - After determining that the paste is completely cured, complete the construction of the core in the usual manner.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

[WARRANTY]

KURARAY MEDICAL INC. will replace any product that is proven to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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"CLEARFIL", "CLEARFIL TRI-S BOND" and "PANAVIA" are trademarks of KURARAY CO., LTD.
"CLEARFIL", "CLEARFIL TRI-S BOND" y "PANAVIA" son marcas comerciales de KURARAY CO., LTD.