
DENTAL ADHESIVE RESIN CEMENT SYSTEM

PANAVIA Veneer LC**ENGLISH INSTRUCTIONS FOR USE****I. INTRODUCTION**

PANAVIA Veneer LC is a light-cure adhesive resin cement system. It consists of the PANAVIA Veneer LC Paste (a light-cure cement paste), CLEARFIL CERAMIC PRIMER PLUS (a silane, zirconia & metal primer), K-ETCHANT Syringe (phosphoric acid), and PANAVIA V5 Try-in Paste. PANAVIA V5 Tooth Primer or CLEARFIL Universal Bond Quick can be chosen for tooth treatment.

The cement paste is a light-cure, resin based material which provides color stability and has radiopacity equal to or greater than 1mm aluminum. It is supplied in Kuraray's ergonomic syringe and dispensed via an angled applicator tip (16G) into an inlay cavity or onto an inlay or a laminate veneer. It is available in 4 shades; Universal (A2), Clear, Brown (A4) and White. It is classified as a Type 2 and Class 2 (Group 1) material by ISO 4049. Shade adaptation may be checked with the PANAVIA V5 Try-in Pastes before cementation.

II. INDICATIONS

Cementation of ceramic and composite inlays, onlays and laminate veneers with less than 2mm thickness.

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers. Patients known to be allergic to any of the ingredients contained in this product.

IV. POSSIBLE SIDE EFFECTS

- The oral mucosa membrane may turn whitish when contacted by the product due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.
- K-ETCHANT Syringe (phosphoric acid) may cause inflammation or erosion due to its chemistry. Use caution to prevent the product from coming in contact with the skin or getting into the eye.

V. INCOMPATIBILITIES

- Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol may cause discoloration and might retard the curing process.
- Do not use hemostatic agents containing ferric compounds since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva due to remaining ferric ions.
- Do not use a hydrogen peroxide solution for cleaning cavities since it may weaken the bonding to the tooth structure.
- Do not use a sodium hypochlorite solution since it may weaken the bonding to the tooth structure.

VI. PRECAUTIONS

- 1. Safety precautions**
 - This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
 - In the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
 - Avoid direct contact with the skin and/or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
 - Exercise caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel or safety glasses to protect them in the event of splashing material.
 - If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>- Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa>- Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol, and rinse with copious amounts of water.
 - Exercise caution to prevent the patient from accidentally swallowing the product.
 - Do not reuse the applicator brush, applicator tip (16G) and needle tip. To prevent cross-contamination, avoid using dispensed liquid/inlay and those brush and tip for different patients. The brush and tip are single-use and should be discarded after using.
- Clean the cavity sufficiently to prevent poor bonding. If the adherent surface is contaminated with saliva or blood, wash it thoroughly and dry. We advise to use KATANA Cleaner for optimal performance of adhesion. Then re-apply the primer or bonding agent.
- This product contains the Ytterbium Trifluoride (less than 15%). Using on children under the age of 6 may have a potential risk of fluorosis.
- Dispose of this product as a medical waste to prevent infection. The tip must be disposed of after covering the tip of the needle to prevent injury.
- The product should be used only for restorations that are less than 2.0 mm thickness and of sufficient translucency. Metal oxide ceramic (e.g. zirconia) should not be used because it may not have sufficient translucency. However, zirconia of our products with high translucency (KATANA Zirconia STML/UTML, KATANA Zirconia Block STML) can be used with a thickness of less than 1.2 mm.

2. Handling and manipulation precautions

- The product must not be used for any purposes other than specified in [II.INDICATIONS].
- The use of this product is restricted to dental professionals.
- Do not use the product as a provisional cement. This material is designed to use as a permanent cement.
- Amalgam or other lining materials remaining in the cavity will prevent the passage of light and the polymerization of the product. Completely remove any lining material when preparing the cavity.
- Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.
- Use a rubber dam to prevent contamination and to control moisture.
- When using temporary materials containing tannin or magnesia, completely remove it to prevent the discoloration.
- When using hemostatic agents containing aluminum chloride, minimize its quantity, and use caution to prevent contact with the adherent surface. Failure to do so might weaken the bond strength to the tooth structure.
- Do not use hemostatic agents containing ferric sulphate. Ferric sulphate may cause the discoloration and weaken the bond strength to the tooth structure.
- Do not mix the product with other materials. The mixed materials may cause a change in physical properties, including a possible decrease in the expected results.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components in the product.
- The liquid bottles should be tightly capped immediately after use to reduce the evaporation of the volatile solvent.
- If the containers and/or instruments for this product are damaged, protect yourself from any danger and immediately discontinue their use.

[PANAVIA Veneer LC Paste]

- Be careful to prevent unnecessary exposure to direct sunlight or operating lights. The paste contains a light-cure catalyst that is highly photo-reactive. During cementation, adjust the angle and/or distance of the dental light to reduce the intensity of light entering the oral cavity to prevent premature polymerization of the paste.
- Excess cement can be removed after tack light-curing it for 1 second. When removing the excess cement, hold the restoration in place to avoid the possibility of lifting the restoration, since there could be some insufficiently cured resin cement.
- When attaching the applicator tip (16G), turn the tip clockwise and attach it securely, thereby preventing paste from leaking at the junction of the applicator tip (16G) and the syringe. In addition, make sure there is no residual paste at the junction of the syringe, which could cause the applicator tip (16G) to fall off.
- Before wiping the residual paste adhering to the applicator tip (16G) or the junction of the syringe with an alcohol gauze pad, squeeze the pad to remove excess alcohol. Use of excessive alcohol in the gauze pad can cause penetration of alcohol into the tip and dilute the paste. In such cases, it may cause a decrease in physical properties from the expected results.
- To expel air from the applicator tip (16G) nozzle and prevent to mix air bubbles to the paste, set the nozzle upward and push the plunger slowly (outside of the patient's mouth) until the paste reaches the nozzle.
- After the paste has been dispensed, the syringe should be capped securely as soon as possible, to prevent the paste from hardening by ambient light, and to prevent foreign matter from entering the syringe.

[PANAVIA V5 Try-in Paste]

- The use of the Try-in Paste should be limited to checking the shade matching with the PANAVIA V5 Paste or PANAVIA Veneer LC Paste.
- The Try-in Paste does not set. Do not use it to cement restorations.
- A transparent liquid may appear at the tip of the Try-in Paste syringe. If this transparent liquid is present, it should be dispensed and discarded, since the separated liquid may affect the shade matching.
- Shade evaluation with Try-in Paste should be performed using the same approximate thickness of Try-in Paste as the hardened cement.
- After use, thoroughly wash the Try-in Paste from the restoration and tooth surface with water to avoid poor adhesion.

[PANAVIA V5 Tooth Primer]

- Use within 5 minutes after dispensing.
- Only use with PANAVIA V5 Paste or PANAVIA Veneer LC Paste.
- If the treated surface is contaminated, wash it with water and dry, or clean with alcohol, or KATANA Cleaner, and treat it again with this primer.

[CLEARFIL CERAMIC PRIMER PLUS]

- Use at once after dispensing. It contains volatile ethanol. As the solvent evaporates, the viscosity increases and may make it difficult to apply.
- Perform cementation soon after treating the restoration surface with this primer.
- If the pre-treated/reated surface is contaminated with saliva or blood, wash it with water, dry, clean with K-ETCHANT Syringe or KATANA Cleaner, and re-treat it with this primer.

[K-ETCHANT Syringe (phosphoric acid)]

- Be careful not to contaminate the treated surface with saliva or blood. If the treated surface is contaminated, re-treat.
- If the product adheres to clothing, wash it off with water.
- After each use, remove the needle tip from the syringe and recap the syringe immediately and tightly.
- When using Tooth Primer, etching vital dentin may cause post-operative sensitivity.

[Dental light-curing unit]

- Use it according to the Instructions for Use for the dental light-curing unit.
- Do not look directly at the light source. Protective glasses are recommended.
- Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible.
- If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Check the conditions required to cure the paste mixture by referring to the light-curing times listed in these Instructions for Use before using the product.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- The product, except for PANAVIA V5 Try-in Paste and K-ETCHANT Syringe, must be stored (2-8°C/36-46°F) when not in use. PANAVIA V5 Try-in Paste and K-ETCHANT Syringe should be stored at (25-27°C/36-77°F). All components stored in the refrigerator should be brought to room temperature for 15 minutes before use in order to restore its normal viscosity and curing performance.
- The product should be kept away from extreme heat, direct sunlight or a flame.
- The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

- Please see the outside of the package for contents and quantity.
- <Principal ingredients>
- PANAVIA Veneer LC Paste: Universal (A2), Clear, Brown (A4), White Silanated spherical silica filler, Urethane dimethacrylate (UDMA), Ytterbium trifluoride, Triethyleneglycol dimethacrylate, Hydrophilic aliphatic dimethacrylate, Hydrophilic amide monomer, Accelerators, di-Camphorquinone, Pigments
 - PANAVIA V5 Try-in Paste: Universal (A2), Clear, Brown (A4), White Glycerol, Silanated colloidal silica, Silanated silica, Colloidal silica, Pigments
 - PANAVIA V5 Tooth Primer 10-Methacryloyloxydicyclo dihydrogen phosphate (MDP), 2-Hydroxyethyl methacrylate, Hydrophilic aliphatic dimethacrylate, Accelerators, Water
 - CLEARFIL CERAMIC PRIMER PLUS 3-Methacryloxypropyl trimethoxysilane, 10-Methacryloyloxydicyclo dihydrogen phosphate (MDP), Ethanol
 - K-ETCHANT Syringe Phosphoric acid, Water, Colloidal Silica, Pigment
 - Accessories Applicator tip (16G)(For Paste), Needle tip (E) (for K-ETCHANT Syringe), Applicator brush (fine <silvers>), Mixing dish

VIII. CLINICAL PROCEDURES

- Cleaning the prepared tooth (tooth, composite resin)** Remove the temporary materials from the prepared tooth in the usual manner, clean the cavity using moisture control. [NOTE] KATANA Cleaner is available to clean the prepared tooth, when contamination occurs. Apply it for 10 seconds with rubbing motion, then rinse with water until the color has been disappeared and dry.
- Trial fitting and adjusting a restoration**
 - Trial fit the prosthetic restoration to check its fit on the prepared tooth, as necessary.
 - If necessary, apply the selected shade of Try-in Paste onto the cementation surface of the restoration and trial fit the restoration onto the prepared tooth. Remove excess Try-in Paste from the margins with a brush. The shades of Try-in Paste correspond to those of the hardened cement paste. Check the shade for best color matching, and then remove the restoration. Using water, completely wash the Try-in Paste from the surface of the restoration and the prepared tooth surface.
- Conditioning the prosthetic restoration surface** Follow the Instructions for Use of the restorative material. In the absence of specific instructions, we recommend the following procedures and application of CLEARFIL CERAMIC PRIMER PLUS:

- If the adherent surface is silica-based ceramic (e.g. conventional porcelain, lithium disilicate), hybrid ceramic or composite resin: Based on the type of restoration, Acid Treatment or Blasting Treatment may be used:
 - Acid Treatment (e.g. conventional porcelain, lithium disilicate): (1) Etch the adherent surface with a hydrofluoric acid solution in accordance with the Instructions for Use of the solution.
 - (2) Rinse the adherent surface with water and dry.

- Blasting Treatment (e.g. hybrid ceramic, composite resin): (1) Roughen the adherent surface by blasting with alumina powder (30-50 µm) using air pressure of 0.1-0.2 MPa (1-2 kgf/cm², 15-29 psi). The air pressure and powder size should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent chipping.
- Clean the prosthetic restoration in an ultrasonic cleaning unit for 2 minutes.

[NOTE]

If the restorations are tried-in after acid treatment or blasting treatment, apply K-ETCHANT Syringe (phosphoric acid) for 5 seconds or KATANA Cleaner for 10 seconds with a rubbing motion then rinse with water until the color has been disappeared and dry the restorations to remove contaminants.

If the adherent surface is metal oxide ceramic (KATANA Zirconia STML/UTML, KATANA Zirconia Block STML):

- Roughen the adherent surface by blasting with alumina powder (30-50 µm) using air pressure of 0.1-0.4 MPa (1-4 kgf/cm², 15-58 psi). The air pressure and powder size should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent chipping.
- Clean the prosthetic restoration in an ultrasonic cleaning unit for 2 minutes followed by drying it with an air stream.

[NOTE]

If the restorations are tried-in after acid treatment or blasting treatment, apply KATANA Cleaner for 10 seconds with a rubbing motion then rinse with water until the color has been disappeared and dry the restorations to remove contaminants.

- Application of CLEARFIL CERAMIC PRIMER PLUS to the prosthetic restoration** Apply the primer to the adherent surface of the restoration with an applicator brush. After application, dry the entire adherent surface sufficiently using mild oil free air flow.

5. Treatment of the prepared tooth

- Choose either PANAVIA V5 Tooth Primer or the CLEARFIL Universal Bond Quick. KATANA Cleaner is available to clean the prepared tooth, when contamination occurs. Apply it for 10 seconds with rubbing motion, then rinse with water until the color has been disappeared and dry.

For PANAVIA V5 Tooth Primer

- If the adherent surface is uncut enamel or when cementing laminate veneers, apply phosphoric acid (K-ETCHANT Syringe) to the enamel surface only and leave it for 10 seconds; rinse and dry.
- Apply Tooth Primer to the entire prepared tooth surface, with an applicator brush and leave it in place for 20 seconds. Use caution not to allow saliva or exudates to contact the treated surfaces.
- Thoroughly dry the entire adherent surface sufficiently by blowing mild, oil-free air. Use a vacuum aspirator to prevent the Tooth Primer liquid from dispersing.

For CLEARFIL Universal Bond Quick

- Choose one of the three etching procedures before applying the bond. If the adherent surface includes uncut enamel or when cementing laminate veneers, application of phosphoric acid (K-ETCHANT Syringe) to the enamel surface is required< (1)-b or (1)-c >.
 - Self-etching procedure
 - Move to below (2) without etching with K-ETCHANT Syringe.
 - Apply K-ETCHANT Syringe to the uncut and/or cut enamel. Leave it in place for 10 seconds, then rinse and dry.
 - Apply K-ETCHANT Syringe to the entire cavity (enamel and dentin), leave it in place for 10 seconds, then rinse and dry.
- Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.
- Apply BOND with a rubbing motion to the entire preparation with the applicator brush. No waiting time is required. Use caution not to allow saliva or exudate to contact the treated surfaces.
- Dry the entire cavity wall sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.
- When cementing laminate veneers, there is no need to light-cure BOND prior to seating the laminate veneers. It should be cured after seating the laminate veneers with the PANAVIA Veneer LC Paste. When cementing inlays or onlays, BOND should be light-cured before seating the restoration. Please confirm the curing time by referencing the following table 1:

Type of light source (Light intensity)	Curing time
High-intensity BLUE LED* (More than 1500 mW/cm ²)	5 sec.
BLUE LED* (800-1400 mW/cm ²)	10 sec.
Halogen lamp (More than 400 mW/cm ²)	10 sec.

*The effective wavelength range of each dental curing unit must be 400 - 515 nm. *Peak of emission spectrum: 450 - 480 nm.

Table 1: Light-curing time of each light source

Restoration type and shade	Maximum restoration thickness	Type of curing light (Light intensity)
Laminate veneer	2.0 mm	High-intensity BLUE LED (≥ 1500 mW/cm ²)
		BLUE LED (1000-1400 mW/cm ²)
Silica-based ceramic (e.g. conventional porcelain, lithium disilicate)	2.0 mm	High-intensity BLUE LED (≥ 1500 mW/cm ²)
		BLUE LED (1000-1400 mW/cm ²)
Hybrid ceramic, composite resin	1.2 mm	High-intensity BLUE LED (≥ 1500 mW/cm ²)
		BLUE LED (1000-1400 mW/cm ²)
Inlays, onlays	2.0 mm	High-intensity BLUE LED (≥ 1500 mW/cm ²)
		BLUE LED (1000-1400 mW/cm ²)
Silica-based ceramic, hybrid ceramics, composite resin	1.2 mm	High-intensity BLUE LED (≥ 1500 mW/cm ²)
		BLUE LED (1000-1400 mW/cm ²)

Table 2: Light-curing time of each Restoration type

Restoration type and shade	Maximum restoration thickness	Type of curing light (Light intensity)		
		High-intensity BLUE LED (≥ 1500 mW/cm ²)	BLUE LED (1000-1400 mW/cm ²)	Halogen lamp (≥ 400 mW/cm ²)
Curing time				
Laminate veneer	2.0 mm	Three times for 3 sec. or Twice for 5 sec.	10 sec.	20 sec.
		High-intensity BLUE LED (≥ 1500 mW/cm ²)	10 sec.	20 sec.
Silica-based ceramic (e.g. conventional porcelain, lithium disilicate)	2.0 mm	Three times for 3 sec. or Twice for 5 sec.	10 sec.	20 sec.
		High-intensity BLUE LED (≥ 1500 mW/cm ²)	10 sec.	20 sec.
Hybrid ceramic, composite resin	1.2 mm	Three times for 3 sec. or Twice for 5 sec.	10 sec.	20 sec.
		High-intensity BLUE LED (≥ 1500 mW/cm ²)	10 sec.	20 sec.
Inlays, onlays	2.0 mm	Three times for 3 sec. or Twice for 5 sec.	10 sec.	20 sec.
		High-intensity BLUE LED (≥ 1500 mW/cm ²)	10 sec.	20 sec.
Silica-based ceramic, hybrid ceramics, composite resin	1.2 mm	Three times for 3 sec. or Twice for 5 sec.	10 sec.	20 sec.
		High-intensity BLUE LED (≥ 1500 mW/cm ²)	10 sec.	20 sec.
Shade and Depth of cure		Universal (A2), Clear, Brown (A4)	1.5 mm	1.0 mm
		White	1.0 mm	

6. Syringe preparation of the cement paste

Select the suitable shade and remove the cap from the selected syringe and attach an applicator tip (16G) securely. Cover the entire syringe with a disposable barrier (e.g. a poly bag) to prevent saliva and blood contamination. Disinfect the syringe by wiping with an absorbent cotton with alcohol both before and after use.

7. Cementing the prosthetic restoration

- Apply the paste over the entire adherent surface of either the prosthetic restoration, or the prepared tooth.
- Place the prosthetic restoration on the prepared tooth.

8. Removing the excess cement

Remove any excess cement using either of the following two methods. Hold the restoration in place during excess removal.
Tack-Cure Clean-Up Technique: Light-cure any excess cement in several spots for 1 second at a distance of 10-15 mm. Hold the prosthetic restoration in position, remove the tack-cured excess cement using a dental instrument. It is advisable to determine in advance the light-curing time of the excess cement by light-curing some paste on a mixing pad.

Wet Clean-Up Technique:

Any excess cement remaining at the margins can be removed with a small brush, foam pellet, dental floss or dental explorer. [For curing margins of the restoration, it is beneficial to cover the margins with a protective gel (e.g. PANAVIA F 2.0 OXYGUARD II) to prevent the formation of an oxygen inhibition layer according to the Instructions for Use.]

9. Final curing

Light-cure the entire surface and margins of the prosthetic restoration that are less than 2.0 mm thick and have sufficient translucency. If the area you want to light-cure is wider than the light emitting tip diameter, extend the exposure process into a few applications. Metal oxide ceramic (e.g. zirconia) should not be used because it may not have sufficient translucency. However, zirconia of our products with high translucency (KATANA Zirconia STML/UTML, KATANA Zirconia Block STML) can be used with a thickness of less than 1.2 mm. Please confirm the curing time by referencing the following table 2:

10. Polishing the margins

Polish the margins using appropriate instruments for the polishing of composite resins.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use the products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

[NOTE]

"PANAVIA", "CLEARFIL" and "OXYGUARD" are trademarks of KURARAY CO., LTD. "KATANA" is a trademark of NORITAKE CO., LIMITED.

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