




KITS, BUR BLOCKS AND ACCESSORIES

PROCEDURAL SOLUTIONS KITS

Tired of Always Searching for the Right Instrument?

Always have the correct instrument-to-task on hand with the SS White® Procedural Solution Kits.

Complete Composite Finishing & Polishing Kit	Complete Restoration Removal Kit	Surgical Procedures Kit	Minimally Invasive Preparation, Finishing & Polishing Kit
Kit Reorder No. 18220 • \$00.00	Kit Reorder No. 18180 • \$00.00	Kit Reorder No. 18170 • \$00.00	Kit Reorder No. 18130 • \$00.00
			
<p>Specifically designed to trim, contour and pre-polish anterior and posterior composite resin materials and polish in one single step.</p> <p>Safe End Trimming & Finishing Burs (1 of each): SE3-10, SE3-20, SE6-10, SE6-20, SE9-10, SE9-20</p> <p>Trimming & Finishing Burs (1 of each): FG 7379, FG 8379, FG 7274, FG 8274</p> <p>Jazz® Polishers: Jazz Supreme Cup, Jazz Supreme Large Flame</p>	<p>Specifically designed for the efficient and atraumatic removal of failed amalgams, composites, crowns and for creating proper endodontic access.</p> <p>Great White® #1 Great White® #2 Revelation® Diamond 806-014C Revelation® Diamond 801-018C Revelation® Diamond 856-016C Great White® Z 801-018 Great White® Z 856-018</p>	<p>This special set of surgical length carbide burs simplifies common surgical procedures including sectioning teeth for extractions, cutting bone and creating proper endodontic access.</p> <p>SL FG #557 SL FG #6 SL FG #8 SL FG #701 SL FG #1702 SL FG #1703 Great White® #2</p>	<p>This kit contains the three different Fissurotomy Burs, which aid in the identification and treatment of hidden caries. Also includes an excellent selection of 12 and 20 blade Finishing Burs including the new anatomically designed FOA1 and the Jazz® Supreme One-Step Composite Polishers.</p> <p>Fissurotomy® Original Fissurotomy® NTF (Narrow Tapered Fissure) Fissurotomy® STF (Shallow Tapered Fissure) FOA 1 FG #7404, FG #8404 Jazz® Polishers: Jazz Supreme Cup, Jazz Supreme Small Flame</p>

UNIVERSAL BUR BLOCKS AND ACCESSORIES

Increase Efficiency and Organization

- The wide color selection allows clinicians to organize their rotary instruments by procedure.
- SS White® Universal Bur Blocks are autoclavable and can be used with diamonds, carbides and polishers.
- The use of bur blocks minimizes damage to instruments and reduces the likelihood of injury from sharp instruments, while creating ideal procedural organization.

SS White® Universal Bur Blocks						
						
Item#	16410 – Blue	16411 – Bordeaux	16412 – Violet	16413 – Green	16414 – Bronze	16415 – Black
PRICE	\$00.00	\$00.00	\$00.00	\$00.00	\$00.00	\$00.00
Each Universal Bur Block holds the following instruments: Front Row: 10 FG/FGSS/SLFG Instruments; Middle Row: 10 FG/FGSS/SLFG Instruments; Back Row: 10 RA/SLRA Instruments						

Autoclavable Bur Blocks & Brushes			
			
Item#	26040 – Wire Cleaning Brush	16392 – FG Thermo Block	16392 – HP/RA Thermo Block
PRICE	\$00.00	\$00.00	\$00.00