

ICE LIGHT

LIGHT TRANSMITTING POST



ICELIGHT is the next generation of aesthetic posts from Danville! **ICE**LIGHT features 70% fill of Zirconia enriched, parallel glass fibers for exceptional strength. Post design is parallel with a tapered end. Its flexibility is close to dentin so stress is distributed, minimizing the risk of cervical fractures.

ICELIGHT is light transmitting for use with light-cured as well as self/dual-cured composites and adhesives. Posts are also presilanated, making placement rapid. The length is easily trimmed eliminating the stocking of multiple lengths. One reamer, per size, gives the ideal post clearance and allows excess composite to escape. Further, there is no composite post with higher radiopacity!

ICELIGHT kit feature 20 posts in sizes 1mm, 1.2, 1.4 & 1.6mm with reamers.

Non-light transmitting, **ICE**POST, are also available. Both are value priced.



ICEPOST Intro Kit 20/pk **90079-01**
 Econo Kit 20/pk 90084-01
 1.0 mm 10/pk Refill 90082-01
 1.2 mm 10/pk Refill 90083-01
 1.4 mm 10/pk Refill 90088-01
 1.6 mm 10/pk Refill 90089-01

ICELIGHT Intro Kit 20/pk **90080-01**
 1.0 mm 10/pk Refill 90378-01
 1.2 mm 10/pk Refill 90380-01
 1.4 mm 10/pk Refill 90382-01
 1.6 mm 10/pk Refill 90384-01

800.827.7940
 TEL (925) 973-0710 FAX (925) 973-0764

Danville
 ENGINEERING & MATERIALS

