

ASTM LEVEL 2 Mask



Skin
Friendly



No Ink



No
Perfume

Features/Benefits:

- ASTM Level 2: Moderate fluid resistance, meant for procedures with moderate fluid/spray
- 3-Ply: Superb breathability without compromising safety
- Pillow Soft Ear loops: Designed for lasting comfort through the workday – no ear fatigue
- No Ink/Perfume/Foam

Size: One Size Fits All
Color: Blue, Green, Pink, White, Yellow
Quantity/UOM: 50pcs/bx, 12bx/cs

Product Code	Description
3100B	ASTM Level 2 Mask, Blue
3100G	ASTM Level 2 Mask, Green
3100P	ASTM Level 2 Mask, Pink
3100W	ASTM Level 2 Mask, White
3100Y	ASTM Level 2 Mask, Yellow

