

ASTM LEVEL 3 Mask



Features/Benefits:

- ASTM Level 3: Used in high fluid/spray procedures. Highest level of protection against fluids and aerosols
- 4-Ply: For extra filtration and maximum fluid protection
- Pillow Soft Ear loops: Designed for lasting comfort through the workday – no ear fatigue
- No Ink/Perfume/Foam

Size: One Size Fits All
Color: Blue, Green, Pink, White, Yellow
Quantity/UOM: 50pcs/bx, 12bx/cs



Product Code	Description
3135B	ASTM Level 3 Mask, Blue

Distributed by: