

COOL BREATHE

PROCEDURE MASK



Level 3



Comfort
Arch



Contour
Fit



SilkSoft

Features

Zero Mask-Fatigue Design

Our Comfort Arch feature paired with our symmetric 4 fold design creates a unique cone shape that ensures extra breathing volume all day long.

Contour Fit

Our vertical contour lines conform to the shape of your face for a secure seal.

SilkSoft Fabric and Ear-Loops

Our masks are hypoallergenic and gentle on the skin and ears to prevent irritation.

ASTM Level 3 Protection

Provides the highest level of fluid protection to keep you and your patients safe.

Size: One Size Fits All

Color: Blue, Black*

Quantity/UOM: 50pcs/bx, 12bx/cs

*Coming Soon

Product Code Description

3145BL

Comfort Arch Mask, Blue

3145BK*

Comfort Arch Mask, Black

