

**INTRODUCING THE ALL NEW!**

**Unlimited** cure depth, no microleakage,

**bulkeZ PLUS™**

does even more!



Unlimited depth of cure



No additional primers or activators needed



Prevents gaps & leakage



One step, no layering or capping saves time



Same great benefits as Bulk EZ, PLUS:



Enhanced shade-matching, polishability, and wear resistance



Do more! Crown fabrication, core build-ups, post cementation



Perfectly balanced flow, stackability, and easy finishing

## Introducing **bulkeZ PLUS™**

The game-changing dual-cure bulk fill composite just got even better, letting you do even more in less time!

Utilizing the same Intellitek Technology for unlimited cure depth and microleakage prevention as the original, Bulk EZ PLUS now incorporates state of the art 80 nm spherical zirconia fillers, enhancing the esthetics with better shade-matching and polishability, and ensuring greater wear resistance. Handling has also been fine tuned to allow for a wider variety of cases, allowing you to do more in less time. Bulk EZ PLUS offers results unlike any other composite from top to bottom!




Learn more about the bulk fill composite that's in a class of its own by calling 800.262.2310 or visiting [Zestdent.com/bulkezplus!](http://Zestdent.com/bulkezplus!)

**ZD ZEST DENTAL SOLUTIONS™**


# THE NEW AND IMPROVED **bulkEZ PLUS™**

- Bulk EZ PLUS™ is an easy-to-place, dual-cure composite that combines unlimited cure depth with flowable cavity adaptation, high strength, and enhanced durability—all in one, simple step.
- Self-cure IntelliTek Technology controls and directs shrinkage, eliminating leakage in all posterior restorations.
- New 80nm spherical zirconia filler particles give market leading esthetics, polishability, and wear resistance.
- Handling and thixotropic properties have been precisely tuned to allow Bulk EZ PLUS to be used in a wider variety of clinical situations.




Images courtesy of Dr. Troy Schmedding, Walnut Creek, CA

## Bulk EZ PLUS is in a Class of It's Own - Here are the Facts:




Bulk EZ PLUS reaches a higher degree of conversion than leading light-cured brands\*



Bulk EZ PLUS offers greater glossiness, polish retention, and wear resistance than leading brands\*

### Bulk EZ PLUS vs Traditional Composite



Light-cured composites shrink towards the light, pulling away and creating leakage. Bulk EZ PLUS eliminates microleakage.\*

## Ordering Information

PART NUMBER	DESCRIPTION
90-00171	Bulk EZ PLUS A1/B1 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00172	Bulk EZ PLUS A2/B2 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00173	Bulk EZ PLUS A3/A3.5/B3 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00174	Bulk EZ PLUS Bleach White Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00175	Bulk EZ PLUS C2/C3 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00176	Bulk EZ PLUS Core White Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00177	Bulk EZ PLUS Trial Kit, 1x 6gm syringe (A1/B1), 1x 6gm syringe (A2/B2), 1x 6gm syringe (A3/A3.5/B3), 1x 6gm syringe (C2/C3), 1x 6gm syringe (BL/OP), 1x 6gm syringe (CORE White), 18x Tips (17 gauge), 18x Tips (19 gauge)
94792	Bulk EZ Mixing Tips (17 Gauge), 25 pack
94793	Bulk EZ Mixing Tips (19 Gauge), 25 pack

Please contact your preferred distributor or Zest representative for additional information on ordering and availability.



**ZEST DENTAL SOLUTIONS**  
2875 LOKER AVENUE EAST  
CARLSBAD, CA 92010  
USA

TEL: 800.262.2310  
FAX: 760.743.7975  
EMAIL: SALES@ZESTDENT.COM  
**WWW.ZESTDENT.COM**