

Great White® Z Diamonds –

**The Hard And Fast Solution For
Cutting Through Zirconia!**



Introducing Great White® Z Diamonds, a superior line of engineered diamond instruments specifically designed to take on the toughest zirconia crowns and copings with ease.

Zirconia crowns and copings are among today's most advanced esthetic ceramics and are gaining in popularity. While they are very strong and esthetically pleasing, they are also extremely difficult and time-consuming to access for endodontic treatment. And, when replacement is necessary, they can be frustrating to remove. Just trimming or adjusting sintered zirconia with standard diamonds and carbide burs can be challenging.

That's why SS White engineered Great White Z Diamonds. Great White Z uses cutting-edge diamond technology to cut easily through this super-hard material. Standard diamonds can only handle the task with great difficulty.



Great White® Z Cuts Zirconia Down To Size

You've always trusted Great White burs for your tough cutting procedures. Now you can be prepared for the day you discover a zirconia coping underneath your patient's porcelain crown...with the diamond advanced enough to carry the Great White name.

Great White Z Features:

- Advanced engineered diamond technology
- Maximum cutting efficiency for zirconia or ceramic crowns, bridges or substructures compared to standard rotary cutting instrumentation (standard diamonds and carbide burs)
- Specialty shapes for slot crown removal technique, endodontic access and adjustment procedures
- Friction grip shank configuration

The Great White Z Diamond Kit contains:

- 2 sizes of round diamonds for use in endodontic access
- 1 round end taper diamond for slot crown removal technique
- 1 football-shaped diamond for adjustment procedures
- Convenient autoclavable bur block







Ordering Information:

Great White Z Introductory Kit

SS White Item #	18160			
Includes 1 each of the following shapes:	GWZ 856-018	GWZ 801-014	GWZ 801-018	GWZ 379-023

Great White Z Refills: 5-Pack

				
SS White Item #	18161	18162	18163	18164
Shape	GWZ 856-018	GWZ 801-014	GWZ 801-018	GWZ 379-023
Description	Round End Taper	Round	Round	Football
Head Length (mm)	8.00	1.30	1.70	4.50
Max Head Diameter (mm)	1.80	1.40	1.80	2.30
LOA (Overall Length mm)	22.00	19.00	19.00	19.00

For optimal results, ensure adequate water flow.

See for yourself how the Great White Z can cut Zirconia down to size. Order your introductory kit today and be prepared for your Zirconia encounter.

To order, contact your dental distributor or call SS White Burs, Inc. at 800-535-2877 for more information.

