

Introducing!



## Sleeve-It™

### Disposable Barrier Sleeves

#### *Help Prevent Cross-Contamination*

Premier's Sleeve-It™ disposable barrier sleeves provide proven protection against cross-contamination. Made from high-grade durable material, Sleeve-It barriers have passed extensive testing<sup>1</sup> to ensure patient safety.

Easy to use, Sleeve-It barriers slide over syringe dispensers and feature a pre-cut opening to pass through the dispensing tip.

Ideal for multiple-use dispensers that cannot be reprocessed using heat-sterilized or high-level disinfectants, Sleeve-It barriers can be conveniently used for each new patient per FDA, CDC and OSHA recommendations.<sup>2\*</sup>



## Sleeve-It™

### Disposable Barrier Sleeves

- ✓ Protection against cross-contamination
- ✓ Passed extensive testing<sup>1</sup> to ensure patient safety
- ✓ Easy to use



**100** SERIES

**9011501** 100 series  
(Fits 0.7 -3.0ml syringe)  
Contains: 300 barrier sleeves



**200** SERIES

**9011503** 200 series  
(Fits 5ml dual-barrel syringe)  
Contains: 300 barrier sleeves



## Vac Attak™ GREEN

### Safe, Effective and EPA-Compliant Evacuation System Cleaner

- ✓ Powerful enzymatic, EPA-compliant cleaner
- ✓ Compatible with amalgam separators
  - pH neutral
  - No chlorine
  - No oxidizing agents
  - Non-foaming
  - Non-corrosive
- ✓ Up to 80 cleanings per jar
- ✓ Economical and easy-to-use powder
- ✓ Fresh citrus scent

**9011105** Six, 800g jars with scoop dispenser



## Brite Shield™

### Instrument Cleaner and Protectant

- ✓ Static soak and ultrasonic cleaner
- ✓ Removes glass ionomer, carboxylate, zinc phosphate cements and bio-burden
- ✓ Unique passivation process coats instruments restricting ruts and corrosion
- ✓ Removes alginate and stone from impression trays
- ✓ Economical and easy to use.  
One jar makes up to 34 gallons

**9011100** Six, 800g jars with scoop dispenser



**premier®** Inspired solutions  
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

<sup>1</sup> Viral penetration: ASTM F 1671  
Synthetic blood penetration: ASTM F 1670  
Resistance to puncture: ASTM F 1342  
Tear resistance: ASTM D 1004  
Tensile testing: ASTM D 882  
Biocompatibility: ISO 10993-5, ISO 10993-10, ISO 10993-12

2. Data on file.

\* IPAC Checklist for Dental Practices-Core Elements, 1st Revision: July 2019 follows FDA recommendation regarding use of a new barrier before each patient on multiple-use dental dispenser devices that cannot be sterilized.

©2020 Premier Dental Co. | All Rights Reserved

8610427 0120064 Rev2 NKP10M