



kuraray

CLEARFIL MAJESTY™ Flow


A Highly Filled, Radiopaque, Micro-Hybrid Composite

A Universal Flowable / And Restorative Solution

- **High radiopacity** for better radiographic detection
- **New syringe with on-demand delivery and bubble-free dispensing**
- **Very high filler loading (81 wt%)** for exceptional strength and wear resistance, even for posterior restorations
- **Liner material or restorative material** with dentin-like properties
- **Lowest polymerization shrinkage** of the current flowables
- **Medium viscosity** so it stays where you place it, yet it flows where it needs to
- **Esthetic shades** for optimal cosmetic results



Precise flow and release formula.



2009
REALITY
Four Star Award
★★★★

THE DENTAL
ADVISOR
++++½


KURARAY AMERICA, INC.

600 Lexington Avenue, New York, NY 10022 | 800.879.1676 | kuraraydental.com | ©2009 Kuraray America, Inc.


Clinical Applications

- Intraoral repairs
- Posterior restorations including Class I's, Class II's and Class V's
- Cavity liner and base for direct restorations, and indirect restorations using the IDS (Immediate Dentin Sealing) technique; the IDS procedure refers to sealing dentin immediately after the preparation and before the final impression

High Radiopacity For Accurate Detection



CLEARFIL MAJESTY™ Flow
290% of aluminum (ISO 4049)




Competitor's Product

Restoration can be done using only CLEARFIL MAJESTY™ Flow



CLEARFIL MAJESTY™ Flow



Competitor's product

CLEARFIL MAJESTY™ Flow is the perfect partner with all restorative materials. It dispenses with virtually no voids or bubbles (other products have voids), and can even be used alone as a final restoration.

Warranty: All data, other than from quoted articles herein, are based on actual measurements performed by Kuraray Medical, Inc. Kuraray Medical Inc. will replace any product that is proved to be defective. Neither Kuraray Medical, Inc. nor the distributor accepts liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user should determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therein. CLEARFIL MAJESTY™ and Kuraray products are available only through Kuraray Authorized dealers. CLEARFIL MAJESTY™ is a trademark of Kuraray Co., LTD. © 2009 Kuraray America, Inc.

CLEARFIL MAJESTY™ Flow Products

Product Name	Item Number	Contents
A1	#2610KA	1 Syringe For Each Shade (3.2g) 15 Needle Tips
A2	#2611KA	
A3	#2612KA	
A3.5	#2613KA	
B1 <i>NEW</i>	#2615KA	
B2 <i>NEW</i>	#2616KA	
C3 <i>NEW</i>	#2618KA	
D3 <i>NEW</i>	#2619KA	
E2	#2620KA	
Cv	#2621KA	
OA3	#2632KA	
BL2	#2635KA	
Needle Tips	#1015KA	20
Shade Guide <i>NEW</i>	#2153KA	1