

**CURING? CANCER SCREENING?  
BEST OF BOTH and MORE...  
FUSION<sup>5</sup> Curing Light**



**Amazing Light for Everyday Practice**

FUSION <sup>5</sup> Bluephase Valo		<p><b>DOE</b> SAVE LIVES!</p> <p><b>Early Detection</b> including HPV</p>
<p><b>Focused</b></p>	<p><b>Dispersed</b></p>	

- 4000+ mW/cm<sup>2</sup>
- 3-sec cure, 10mm depth
- Super efficient and low heat
- Caries, canal and cancer detection
- Cold laser therapy

**You've Got the Power and More**



FUSION Curing Light is a 21st century light that is versatile, cost-effective, and valuable in everyday practice. Adding to the value is the ability to conduct oral cancer examination.

**Dr. Howard Glazer**



The Fusion Light is one GREAT all around light for the dental office. You get the benefits of three lights in one! Curing, Transillumination and Oral Cancer Screening. The Curing power is remarkable.

**Dr. John Comisi**




**Version 5**

**DentLight**

# FUSION<sup>5</sup>

## Amazing Power and Style in the Palm of Your Hand!

- Unique focused beam and advanced backlit LCD display with smart digital programs:
  - 4000 mW/cm<sup>2</sup> Plasma Power mode for fast/deep cure
  - Pulsed 2000 mW/cm<sup>2</sup> mode for low heat/shrinkage cure
  - Built-in oral exam, treatment and bleaching modes
  - Battery and thermal bars for smart control
- Wide spectrum of 415 - 490nm to cure all materials
  - Options for lower wavelengths down to 385nm
- Sleek metal keypad with large activation and program buttons
- Modular design with easy upgrade and service parts
- Cordless with direct power plug that eliminates downtime
- 90 minutes fast recharge time
- Super slim and ergonomic
  - Tip height only 11.5mm
  - Weighs only 3.9oz
  - Measures at 7.8 inch




### Built to Last

Robust anodized aluminum metal construction distributes low LED heat, prevent scratches, discoloration and damage from accidental falls.

### Lowest Profile


360° swivel head for easier maneuvering and access




### Smart Lithium-Ion Battery You Can Count On!

- Cordless power with 70 minutes plasma power run time
  - 1400 3-second cures
- Multiple charge options
  - Space saving counter stand
  - Direct plug in
- Constant output throughout battery run time

### Curing Light Dispersion



### ADD-ON FUNCTIONS



Cancer Screening




Adhesive Views



Caries Detection



Photodynamic Therapy/Disinfection



Chairside Bleaching

### 7800080 FUSION 5 Curing Light Kit:

Including Handpiece, Charging Stand, Curing Caps (tacking, 9mm and 13mm), Power Adapter, Light Shield, Barrier Sleeves (100), and Instructions for Use.

**2-year warranty** (1-year Battery and Charging Stand)

**DentLight**

**CALL 800-763-6901 today**

### DENTLIGHT INC

1825 Summit Ave. #210

Plano, TX 75074, US

Fax: 972-346-6550

Email: [sales@dentlight.com](mailto:sales@dentlight.com)