


JM Disposable Dental Needles


- Indicator dot on hub provides easy identification of the lancet bevel position for less insertion pain at the injection site
- Lancet bevel point at cartridge end prevents blockage of anesthetic and ensures easy flow of anesthetic solution
- Universal plastic hub fits most syringes
- Color coding for easy identification
- Economical



Lancet bevel prevents rubber chips from blocking needle bores

One dot needle

Designed for both Luer-Lock and screw-on type syringes



Ordering Information

(100) sterile, disposable dental needles are included in each box.

| | |
|--------------------------|---------|
| 25 G Long (30 mm) | 20-25GL |
| 27 G Short (21 mm) | 20-27GS |
| 27 G Long (30 mm) | 20-27GL |
| 30 G Short (21 mm) | 20-30GS |
| 30 G Long (25 mm) | 20-30GL |

