


Prelude^{one}™

One Component Universal Adhesive

Prelude One is a light-curable single component self-etch adhesive designed to bond composite resins to most dental surfaces. These surfaces include dentin, enamel, zirconia, alumina, metals, porcelain, lithium disilicate, and composite restorations. As a universal adhesive, Prelude One is truly compatible with light cured, dual cured, and auto cured composites (no activator required).

Shear Bond Strength on Tooth Surface (LC Mode)




■ SBS to cut enamel
■ LC (MPa)

■ SBS to dentin
■ LC (MPa)



Shear Bond Strength on Other Surface (LC Mode)



■ SBS to Lithium disilicate (E-Max CAD, Ivoclar) with HF treatment (MPa)

■ SBS to precious metal (Panceram 62, Pentron) MPa

■ SBS to non-precious metal (Cocr) MPa


■ SBS to Zirconia (3M Lava) MPa

Available in 5ml bottle and Unitdose

ADVANTAGES

- **Convenient one step application**
- **For use on all surfaces**
- **No post-op sensitivity**
- **Extremely high bond strengths to all direct and indirect restorations**
- **Compatible with dual and self cure composites without the use of an additional activator**
- **5 1/2 micron film thickness**

Shear Bond Strength on Tooth Surface (AC Mode)



For more information visit our website or call our Customer Service Hotline.

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

94196 Rev A