



# TwinPower Turbine® Basic

Air Driven, High Speed Handpiece




- High power, up to 22 watts
- Safe, rapid braking
- Chrome body
- Ceramic bearings
- 4-hole connection



# TwinPower Turbine Basic

The new TwinPower Turbine Basic is an excellent choice for any practice. It offers several features of the original product line, above and beyond those in its class.

Power	Safety	Design
The Basic model offers the same cutting efficiency as the TwinPower Turbine standard head. Its double-impeller rotor design produces up to 22 watts of power.	A pressurized air system prevents suck back in the air line resulting in superior infection control. Additionally, a unique rubber brake ring design stops the turbine within 2 seconds, allowing for safer and more efficient preparations.	The body of the handpiece is chrome and comes with a convenient push-button chuck. Lightweight, at just 45 grams, it also helps reduce operator fatigue.

Ordering Information and Specifications			
Specifications	 Ceramic ball bearing handpiece	 Push-button chuck	 4 Hole
Name	TwinPower Turbine Basic		
Model	PAR-4HEX-B	Body	Chrome
Product number	16-5360838	Zero suck back	Yes
Kit includes	Lubricant spray nozzle	Quick brake	Yes, within 2 seconds
Power (watts)	22	Water check valve	Yes
Rotation speed (rpm)	370,000 / ± 30,000	Dental washer disinfectant	Yes
Water spray ports	3	Air flow	37±3NL/min
Head diameter (mm)	10.5	Driving air pressure	0.2 – 0.29 MPa/29 - 42 psi
Head height (mm)	13.2	Optics	No
Weight (g)	45	Warranty	1 year